



# PRODUCTS FOR SALIVARY DUCT ACCESS

## Salivary Duct Access

Gain and maintain access to the salivary duct while reducing trauma and maximizing efficiency.

Sialendoscopy represents an exciting **minimally invasive** alternative to a standard open surgical procedure to treat salivary duct disorders. Sialendoscopy has a variety of advantages for the patient, including minimal recovery time, no incision or subsequent scarring, minimized potential for nerve damage, preservation of the salivary gland<sup>1</sup> and early return to a normal diet. Sialendoscopy allows **stones to be removed in the office** with a low risk of major complications to the patient.

Cook Medical continually collaborates with physicians to improve procedures. Cook has partnered with leading physicians to create salivary duct access products that are designed to maintain access throughout the procedure and work together for greater control. This **control** can help reduce trauma for the patient and increase efficiency with every step.

“Cook’s salivary duct access tools allow me to gain and maintain access to the duct, which gives me much greater control throughout the entire procedure.”

Jack Kolenda, MD  
Chief, Department of Otolaryngology  
Oakville-Trafalgar Memorial Hospital  
Oakville, Ontario

Dr. Kolenda is the inventor of the Kolenda Salivary Access Introducer Set and receives a royalty payment based upon Cook’s license to use this technology.

### Reference

1. Wallace E, Tauzin M, Hagan JP, et al. Management of giant sialoliths: review of the literature and preliminary experience with interventional sialendoscopy. *Laryngoscope*. 2010;120(10):1974-1978.

### Resources

1. Chossegros C, Guyot L, Richard O, et al. A technical improvement in sialendoscopy to enter the salivary ducts. *Laryngoscope*. 2006;116(5):842-844.
2. Walvekar R, Razfar A, Carrau RL, et al. Sialendoscopy and associated complications: a preliminary experience. *Laryngoscope*. 2008;118(5):776-779.

- AORTIC INTERVENTION
- CRITICAL CARE
- ENDOSCOPY
- INTERVENTIONAL RADIOLOGY
- LEAD MANAGEMENT
- PERIPHERAL INTERVENTION
- SURGERY
- UROLOGY
- WOMEN’S HEALTH

### Customer Service Centers

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)




Illustration by Lisa Clark


# GAIN AND MAINTAIN ACCESS

1

## WIRE GUIDE

Gain and maintain access with a wire guide.

Sialendoscopy starts with accessing the salivary duct through the papilla. A flexible Cook soft-tip wire guide can be easily maneuvered into the duct. This guide allows you to easily gain and maintain ductal access throughout the procedure, reducing trauma and maximizing efficiency.




2

## DILATOR SET

Minimize trauma with sequential dilation.

Once you establish access using the Cook soft-tip wire guide, a series of flexible dilators can be introduced over the wire to easily expand the papilla and prepare the salivary duct for the introduction of procedural instruments.




3

## KOLENDA INTRODUCER SHEATH

Create a true working channel.

After the dilator set is used to expand the papilla, the sheath can be passed over the wire to create an open working channel into the salivary duct. The sheath also serves to protect the ductal wall, allowing easy insertion and removal of the procedural instruments.



Illustrations by Lisa Clark

### WIRE GUIDE

Order Number	Reference Part Number	Diameter inch	Length cm	Tip Configuration
G02255	C-SF-18-40	.018	40	straight

ADAPTOR	G05301	PMLLA-VTA-L
---------	--------	-------------

### SALIVARY ACCESS DILATOR SET

Order Number	Reference Part Number	Dilator OD Fr/mm/inch
G24131	SDDS-18-20-HC	4.0/1.35/.053
		5.0/1.67/.066
		6.0/2.0/.079
		7.0/2.33/.092

### KOLENDA SALIVARY ACCESS INTRODUCER SET

Order Number	Reference Part Number	Introducer ID mm/inch	Wire Guide Diameter inch	Sheath Length cm
G24132	SDAS-5.0-18	1.9/.074	.018	5.5
G24129	SDAS-6.0-18	2.2/.087	.018	5.5