



# Van Schie

BEACON® TIP SEEKING CATHETER




## Simplify cannulation.

An Accessory for Endovascular Procedures



# Van Schie

## BEACON® TIP SEEKING CATHETER



The Beacon Tip Van Schie Seeking Catheter is used to ease navigation during the catheterization of bifurcated endoluminal AAA stent grafts. Different available tip designs allow the Beacon Tip Van Schie to simplify cannulation within iliacs of varying size and tortuosity.

### Compatibility

The Beacon tips have five different configurations from which to choose, each sold separately. The Van Schie catheter accommodates a .035 inch wire guide and comes standard in 5 Fr with a length of 65 cm.

### Visibility & Durability

The radiopaque Beacon tips enhance visualization during cannulation. Their nylon-braided construction provides torque control, as well as enhanced flexibility and durability.

Global Product Number	Order Number	Fr	Accepts Wire Guide Diameter inch	Length cm
G13795	HNBR5.0-35-65-P-NS-VANSCHIE1	5.0	.035	65
G13794	HNBR5.0-35-65-P-NS-VANSCHIE2	5.0	.035	65
G13792	HNBR5.0-35-65-P-NS-VANSCHIE3	5.0	.035	65
G13783	HNBR5.0-35-65-P-NS-VANSCHIE4	5.0	.035	65
G13793	HNBR5.0-35-65-P-NS-VANSCHIE5	5.0	.035	65



[www.cookmedical.com](http://www.cookmedical.com)

#### COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.  
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

#### COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA  
Phone: 905.640.7110, Toll Free: 800.668.0300

#### WILLIAM A. COOK AUSTRALIA PTY. LTD.

95 Brandl Street, Brisbane Technology Park, Eight Mile Plains  
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

#### WILLIAM COOK EUROPE ApS

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2009

AI-BM-VANDS-EN-200910

AORTIC  
INTERVENTION

CARDIOLOGY

CRITICAL  
CARE

ENDOSCOPY

INTERVENTIONAL  
RADIOLOGY

PERIPHERAL  
INTERVENTION

SURGERY

UROLOGY

WOMEN'S  
HEALTH