

### DEVICE PLANNING AND SIZING WORKSHEET

Date: \_\_\_\_\_ Patient ID: \_\_\_\_\_

Hospital: \_\_\_\_\_

Physician name: \_\_\_\_\_

Physician phone #: \_\_\_\_\_

Physician e-mail: \_\_\_\_\_

Physician signature: \_\_\_\_\_


### Anatomical measurements

Introduction site

☐ Right iliac External iliac (EI) measurement \_\_\_\_\_mm

☐ Left iliac External iliac (EI) measurement \_\_\_\_\_mm

Introduction system: 28-34 mm grafts = 7.7 mm (OD)  
36-42 mm grafts = 8.5 mm (OD)



Diameters: Dp  Dd

Lengths: AL

Aortic vessel diameters; CT diameters are outer wall to outer wall

### Aneurysm length (AL)

If AL is > 30 mm with appropriate proximal and distal neck lengths, see the ordering chart for recommended proximal and distal grafts. If AL is ≤ 30 mm with appropriate proximal and distal neck lengths, use a single proximal graft.

### Device overlapping

A minimum overlap of 2 stents (50 mm) is required; a 3-4 stent (75-100 mm) overlap is recommended. Ensure that the distal graft does not extend beyond the top sealing stent of the proximal graft.

### Tapered anatomy

If Dp (diameter, proximal neck) is 4 mm > Dd (diameter, distal neck), then use a proximal tapered graft.

### Proximal grafts

Proximal aortic vessel diameter (Dp) =

AL plus  25 =  Proximal working length  
Landing zone

Select the graft diameter from the vessel oversizing chart below. Select the graft length that is closest to the working length.

#### Proximal Graft (ZTEG-2P)

| Diameter mm | Length mm          |
|-------------|--------------------|
| 28, 30, 32  | 120, 140, 200      |
| 34, 36, 38  | 127, 152, 202      |
| 40, 42      | 108, 135, 162, 216 |

#### Proximal Tapered Graft (ZTEG-2PT)\*

| Diameter mm | Length mm |
|-------------|-----------|
| 32          | 160, 200  |
| 34          | 157, 197  |
| 36          | 157, 197  |
| 38          | 152, 202  |
| 40          | 158, 208  |
| 42          | 158, 208  |

\*Proximal tapered graft diameter is 4 mm > the distal sealing stent.

☐ Proximal order number = ZTEG-2P -  -  - PF-US  
Graft diameter Graft length

☐ Proximal tapered order number = ZTEG-2PT -  -  - PF-US  
Graft diameter Graft length

### Distal grafts

Distal aortic vessel diameter (Dd) =

AL plus  25 =  Distal working length  
Landing zone

Select the graft diameter from the vessel oversizing chart below. Select the graft length that is closest to the working length.

#### Distal Graft (ZTEG-2D)

| Diameter mm | Length mm |
|-------------|-----------|
| 28, 30, 32  | 147, 207  |
| 34, 36, 38  | 136, 186  |
| 40, 42      | 144, 198  |

Distal order number = ZTEG-2D -  -  - PF-US  
Graft diameter Graft length

### Diameter

Distal diameter of proximal graft must be the same as proximal diameter of distal graft.

### Graft-to-vessel oversizing (Dp and Dd)

Apply oversized graft calculation to the sizing chart when selecting graft diameter.


#### Vessel oversizing chart (10-25%)\*\*

| Aortic Vessel Diameter mm | Graft Diameter** mm |
|---------------------------|---------------------|
| 24                        | 28                  |
| 25, 26, 27                | 30                  |
| 28, 29                    | 32                  |
| 30                        | 34                  |
| 31, 32                    | 36                  |
| 33, 34                    | 38                  |
| 35, 36                    | 40                  |
| 37, 38                    | 42                  |

\*\*Additional considerations may affect the choice of diameter.

Patient ID: \_\_\_\_\_

## CASE NOTES



### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com/web/aortic-intervention](http://www.cookmedical.com/web/aortic-intervention)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

