

Intubate with confidence.



Frova
INTUBATING CATHETER



www.cookmedical.com

Frova

INTUBATING CATHETER

The Frova Intubating Catheter is used to facilitate endotracheal intubation in patients when visualization of the glottis is inadequate.


The Frova Intubating Catheter offers many benefits:

- The blunt, angled tip both assists in proper placement beyond the vocal cords and minimizes trauma to internal structures.
- The removable Rapi-Fit® adapter permits the use of a ventilation device, if necessary.
- The optional stiffening cannula adds rigidity to all but the distal 5 cm of the catheter.
- The through-lumen design and distal sideports ensure adequate airflow.
- The centimeter marks facilitate accurate placement with endotracheal tubes.




FROVA INTUBATING INTRODUCER

Radiopaque




STIFFENING CANNULA

(This cannula is included only with -FII Reference Part Numbers.)



RAPI-FIT ADAPTER
15 mm connector



RAPI-FIT ADAPTER
Luer lock connector

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Catheter ID mm	For Placement of ET Tubes ¹ with ID mm
G48303	C-CAE-14.0-70-FIC ²	14.0	70	3.0	≥ 6
G53517	C-CAE-14.0-70-FIC-SPOPS	14.0	70	3.0	≥ 6
With Rapi-Fit Adapters					
G23645	C-CAE-14.0-70-FI	14.0	70	3.0	≥ 6
With Rapi-Fit Adapters and Stiffening Cannula					
G12591	C-CAE-8.0-35-FII	8.0	35	1.6	≥ 3
G23646	C-CAE-14.0-70-FII	14.0	70	3.0	≥ 6

1. An endotracheal tube is not included.

2. This product is supplied in a single box that contains 10 catheters with individually labeled packages.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details. For complete instructions, contraindications, warnings, and precautions see the Instructions for Use that are included with the product.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 972519996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com

