




ENDOBONCHIAL BLOCKER SETS




Select a simple alternative
for one-lung ventilation.

Cook's endobronchial blockers allow one-lung ventilation through a conventional single-lumen endotracheal tube (SLT), which is easier to place and poses less risk of trauma to internal structures than a double-lumen tube (DLT).¹ Loss of the airway, a potential problem during the exchange of an SLT for a DLT, is also eliminated with this method.²



Arndt

ENDOBONCHIAL BLOCKER



FEATURES & BENEFITS

- Adjustable guide loop coupled with an appropriately sized bronchoscope enables precise placement.
- Low-pressure, high-volume balloon creates excellent surface area contact with the inner bronchial wall while minimizing potential trauma to the bronchus.
- Distance markings show how far the blocker has been advanced.
- Suction adapter included with set offers another method of lung deflation.


Order Number	Reference Part Number	Fr	Length cm
Arndt			
G44109	C-AEBS-5.0-50-SPH-AS	5.0	50
G44114	C-AEBS-7.0-65-SPH-AS	7.0	65
G44120	C-AEBS-9.0-78-SPH-AS	9.0	78

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



Cohen

ENDOBONCHIAL BLOCKER



FEATURES & BENEFITS


- Unique deflecting tip provides wide range of motion to precisely direct the blocker under bronchoscopic guidance.
- Indicator arrow on the turning wheel shows which direction the tip deflects, ensuring reliable access to the desired lung.
- High-traction, silicone grip on introducer shaft increases operator control, resulting in more precise device placement.
- Low-pressure, high-volume balloon creates excellent surface area contact with the inner bronchial wall while minimizing potential trauma to the bronchus.
- Distance markings show how far the blocker has been advanced.
- Suction adapter included with set offers another method of lung deflation.

Order Number	Reference Part Number	Fr	Length cm
Cohen			
G44122	C-AEBTS-9.0-65-SPH-AS	9.0	65

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

REFERENCES

1. Knoll H, Ziegeler S, Schreiber J, et al. Airway injuries after one-lung ventilation: A comparison between double-lumen tube and endobronchial blocker. *Anesthesiology*. 2006;105(3):471-477.
2. Arndt GA, DeLessio ST, Kranner W, et al. One-lung ventilation when intubation is difficult—presentation of a new endobronchial blocker: a randomized, prospective, controlled trial. *Acta Anaesthesiol Scand*. 1999;43(3):356-358.



AORTIC INTERVENTION

CRITICAL CARE

ENDOSCOPY

INTERVENTIONAL RADIOLOGY

LEAD MANAGEMENT

PERIPHERAL INTERVENTION

SURGERY

UROLOGY

WOMEN'S HEALTH



Customer Service Centers

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com
www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,
E-mail: orders@cookmedical.com

Australia: ABN: 51 141259118
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,
E-mail: wcacustserv@cookaust.com.au

