


Expand your ability to treat obstructive salivary gland disorders.



SialoCath™

SALIVARY DUCT CATHETER

Focused

This catheter is designed specifically for use in the salivary ducts.


Versatile

- Use it in the submandibular or parotid glands.
- Use it to irrigate salivary ducts or flush stone fragments.
- Insert it over a wire guide or through a Kolenda Introducer Sheath.

Atraumatic

This flexible, soft-tipped catheter supports minimally invasive intervention and provides a gentle treatment option for in-office and surgical procedures.






SialoCath™

SALIVARY DUCT CATHETER

The SialoCath Salivary Duct Catheter is intended to aid in the removal and treatment of salivary duct stones and other duct conditions through irrigation of the ducts and/or through flushing of the salivary ducts. This flexible, single-lumen catheter can be inserted into the salivary duct either over a wire guide or through an appropriately sized introducer sheath.

Order Number	Reference Part Number	Catheter OD mm/inch	Catheter ID mm/inch	Wire Guide Diameter inch	Catheter Length cm	Tip Configuration
G35166	SDC-2.2-20	.74/.029	.53/.021	.018	15	straight


Related Products



Kolenda Introducer Sheath

Intended for the introduction of diagnostic or interventional devices in the submandibular salivary duct.

Order Number	Reference Part Number	Introducer ID mm/inch	Wire Guide Diameter inch	Sheath Length cm
G24129	SDAS-6.0-18	2.2/.087	.018	5.5
G24132	SDAS-5.0-18	1.9/.074	.018	5.5



Wire Guides

Used to gain and maintain access to the salivary duct.

Order Number	Reference Part Number	Wire Guide Diameter inch	Wire Guide Length cm	Tip Configuration
G02274	C-SF-15-50	.015	50	straight
G02255	C-SF-18-40	.018	40	straight

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com

