


Transfer with ease.



Guardia™ Access Nano

CURVED EMBRYO TRANSFER CATHETER

Atraumatic


The reduced outer diameter of the bulb tip is designed to enable atraumatic catheter passage through the cervix.

Flexible

The soft transfer catheter provides the flexibility to navigate different patient anatomies.



www.cookmedical.com



Guardia Access Nano

Cook Medical now offers another choice for embryo transfer: the Guardia Access Nano. This precurved embryo transfer catheter has been optimized to ease passage through the cervix.

Cook's Guardia Access Nano is the result of refining a best-selling product to suit the comfort and needs of patients. It has Cook's traditional distal curve and bulb tip, which has a reduced outer diameter to address the challenges of smaller anatomy. The Guardia Access Nano's transfer catheter also has an internal support cannula, which is more flexible during embryo placement than an external support cannula. The transfer catheter is available in two lengths, 24 cm and 30 cm.

Order Number	Reference Part Number	Guide Catheter Fr	Guide Catheter Length cm	Transfer Catheter Fr	Transfer Catheter Length cm
G24216	K-JETS-551910-S	5.5	17.3	2.8	24
G24215	K-JETS-551910-L	5.5	17.3	2.8	30

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Guardia Obturator

Use Cook's Guardia Obturator to assist an approved embryo transfer catheter in gaining uterine access to place in vitro fertilized (IVF) embryos into the uterine cavity.

Order Number	Reference Part Number	OD inch	Length cm
G24217	J-UOB-2828-ST	.035	28

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com

