


Achieve visualization of the uterus and fallopian tubes.



SILICONE BALLOON HSG CATHETER

Control

The small diameter of the catheter is designed to allow easy passage through the cervix.


Visibility

The size of the balloon provides an enhanced view during evaluation.

Versatility

The catheter can be used to perform SHG and HSG procedures.





SILICONE BALLOON HSG CATHETER

The Cook® Silicone Balloon HSG Catheter is used to access the uterine cavity for saline infusion sonohysterography or to inject contrast medium into the uterine cavity and fallopian tubes for evaluation of tubal patency.

- Silicone construction; contains no latex.
- Available with or without sideport.
- Includes an appropriately sized syringe.

Cook Medical—Developing simple solutions for reproductive health.

Order Number	Reference Part Number	Fr	Length cm	Balloon Volume mL	Syringe Volume mL
Closed Tip with Sideport					
G17472	J-CHSG-503001	5.0	30	1.0	1.0
G17471	J-CHSG-703001	7.0	30	1.5	3.0
Open Tip					
G17147	J-CHSG-503000	5.0	30	1.0	1.0
G17870	J-CHSG-703000	7.0	30	1.5	3.0

Sold in boxes of 10.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com

