


Selectively evaluate the fallopian tubes under direct vision.



Modified Novy

CORNUAL CANNULATION SET

Clarity

The clear tip provides an enhanced view of the catheter's pathway.

Precision


Direct vision through a hysteroscope enables careful placement.

Uniqueness

The set is one of the only products that is available for hysteroscopy of the proximal fallopian tubes.



www.cookmedical.com




Modified Novy

CORNUAL CANNULATION SET

Intended for use through the operating channel of a hysteroscope or other uterine access device for hysteroscopic or fluoroscopic selective catheterization and cannulation of the proximal fallopian tube(s). This may be followed by the introduction of chromotubation solution or contrast medium to evaluate tubal patency. The wire guide in the Modified Novy Cornual Cannulation Set is intended only to facilitate placement of the inner catheter.¹

- Torcon® polyethylene introducing catheter facilitates access to the cornual os.
- Transparent catheter tip provides maximum visualization of the inner catheter and wire guide.

1. The Novy is not intended for complete catheterization of the fallopian tube(s) or for tubal recanalization. Also, it is not intended for embryo transfer or other assisted reproduction procedures.



Order Number	Reference Part Number	Introducing Catheter Fr	Introducing Catheter Length cm	Inner Catheter Fr	Inner Catheter Length cm	Wire Guide Diameter inch/mm
Straight Introducing Catheter						
G26822	J-NCS-503500	5.0	35	3.0	50	.018/.46
Curved Introducing Catheter						
G17478	J-NCS-503570	5.0	35	3.0	50	.018/.46
G17558	J-NCS-504070	5.0	40	3.0	55	.018/.46

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 97251996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com

