

Maintain access
at all times.



Fusion[®]OASIS[®]

ONE ACTION STENT INTRODUCTION SYSTEM

DUAL PLATFORM ERCP

ACCESS

STONE MANAGEMENT

STRICTURE MANAGEMENT



www.cookmedical.com



Fusion® OASIS®

ONE ACTION STENT INTRODUCTION SYSTEM

Used for endoscopic biliary stent placement to drain obstructed bile ducts.

- Eliminate the need for recannulation when placing multiple stents, repositioning a stent, or fully removing the stent from the ductal system
- Integrated guiding and pushing catheter provides optimal pushability when negotiating tortuous anatomy
- Compatible with a wide variety of our exceptional plastic stents

Order Number	Reference Part Number	Pushing Catheter Fr	Guiding Catheter Fr	Pushing/ Guiding Catheter Length cm	Wire Guide Diameter inch	Minimum Accessory Channel mm
G31526	FS-OA-8.5	8.5	5	170/200	.035	4.2
G31527	FS-OA-10	9	5.5	170/200	.035	4.2

Not for use with pigtail stents. 1 radiopaque band at IDE port.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Fusion Pushing Catheter

Order Number	Reference Part Number	Pushing Catheter Fr	Length cm	Minimum Accessory Channel mm
G31550	FS-PC-5	5	170	3.2
G31551	FS-PC-7	7	170	3.2

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 97251996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com

