


# New possibilities for cannulation.



## Fusion<sup>®</sup> LoopTip<sup>™</sup>

WIRE GUIDE

DUAL PLATFORM ERCP


ACCESS

STONE MANAGEMENT

STRICTURE MANAGEMENT



[www.cookmedical.com](http://www.cookmedical.com)



## Fusion® LoopTip™

### WIRE GUIDE

Used to assist in cannulation of the biliary and pancreatic ducts and to aid in bridging difficult strictures during ERCP.

- Closed distal loop delivers potentially less traumatic access and navigation of the pancreatobiliary ductal system
- Radiopaque coil spring provides excellent fluoroscopic visibility
- Spiral markings extending to the distal end of the wire guide provide endoscopic view of wire guide movement during cannulation
- Nitinol core for kink resistance and shape retention

| Order Number | Reference Part Number | Wire Guide Length<br>cm | Wire Guide Diameter<br>inch | LoopTip Size<br>mm |
|--------------|-----------------------|-------------------------|-----------------------------|--------------------|
| G50975       | FS-LT-35-205          | 205                     | .035                        | 2x4                |
| G50976       | FS-LT-35-260          | 260                     | .035                        | 2x4                |
| G50977       | FS-LT-35-480          | 480                     | .035                        | 2x4                |

May be left in place during sphincterotomy. **Caution:** Not for use with metal tip devices.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)

Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)

Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)

Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)

Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)

Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)

France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)

Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)

Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)

Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)

Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)

Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)

Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)

Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)

Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)

Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)

Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)

Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)

United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

