

# Plastic Stents

BILIARY AND PANCREATIC

Established Clinical Performance



[www.cookmedical.com](http://www.cookmedical.com)

## TABLE OF CONTENTS

### Biliary Stents

|                                                                |    |
|----------------------------------------------------------------|----|
| Zimmon Double Pigtail Stent .....                              | 1  |
| Compass BDS Double Pigtail Stent.....                          | 3  |
| Solus Double Pigtail Stent with Introducer.....                | 5  |
| Cotton-Leung Stent .....                                       | 7  |
| Cotton-Leung Sof-Flex Stent .....                              | 9  |
| Cotton-Huibregtse Stent.....                                   | 11 |
| ST-2 Soehendra Tannenbaum Stent.....                           | 13 |
| OASIS One Action Stent Introduction System -<br>Preloaded..... | 15 |

### Pancreatic Stents


|                                   |    |
|-----------------------------------|----|
| Geenen Stent .....                | 17 |
| Geenen Sof-Flex Stent .....       | 19 |
| Zimmon Single Pigtail Stent ..... | 21 |
| Johlin Wedge Stent .....          | 23 |

### Introducers

|                                                           |    |
|-----------------------------------------------------------|----|
| Fusion OASIS One Action Stent<br>Introduction System..... | 25 |
| OASIS One Action Stent<br>Introduction System.....        | 27 |
| Other Stent Introducers .....                             | 29 |

|                         |    |
|-------------------------|----|
| Order Information ..... | 31 |
|-------------------------|----|

# Biliary Stents




## BILIARY STENTS

# Zimmon<sup>®</sup>

## BILIARY STENTS


When performance is essential, choose the Zimmon Biliary Stent. This double pigtail design minimizes stent migration while side holes enhance duct drainage.



# Zimmon®

### BILIARY STENTS

- Minimizes migration with double pigtail configuration
- Smooth cannulation facilitated by tapered tip
- Pushability and softness delivered by polyethylene material
- Simplify inventory management and procedure efficiency with stent sets, which include pushing catheter




*Image courtesy of Anil Sharma, MD,  
Highland Hospital, Rochester, NY*

## BILIARY STENTS

# Compass BDS™

### BILIARY PLASTIC STENT

With the Compass BDS, you can place biliary stents with ease and accuracy. This bidirectional stent, which can be simply loaded and delivered, provides enhanced drainage for your patients.



# Compass BDS™

### BILIARY PLASTIC STENT

- Ease of loading and introduction into the ductal system is facilitated by the dual tapered ends
- Enhanced drainage with 17 drainage holes for each pigtail
- Fluoroscopic visualization for positioning is enhanced with two radiopaque markers




*Image courtesy of Dr James W. Ostroff,  
UCSF Medical Center, San Francisco, CA*

## BILIARY STENTS

# Solus<sup>®</sup>

### DOUBLE PIGTAIL STENT WITH INTRODUCER


The Solus Double Pigtail Stent is supple yet strong and completely radiopaque, making it an innovative solution for biliary drainage. For the ultimate in flexibility and strength, choose the Solus solution.



### Solus®

#### DOUBLE PIGTAIL STENT WITH INTRODUCER

- Minimizes migration with double pigtail configuration
- Fluoroscopic visibility enhanced by Sof-Flex, a soft pliable radiopaque material
- Accurate position is achieved by visibility of radiopaque markers.
- Efficient drainage provided by multiple side holes
- Simplify inventory management and enhance procedure efficiency—with included one-step introducer




*Image courtesy of Dr. Wahid Wassef, University of  
Massachusetts Medical Center, Worcester, MA*

## BILIARY STENTS

# Cotton-Leung<sup>®</sup>

### BILIARY STENT


The subtle, natural curve of the Cotton-Leung Biliary Stent enhances its introduction and positioning capabilities, making it the stent of choice for a variety of clinical needs.



# Cotton-Leung®

## BILIARY STENT

- Minimizes migration with proximal and distal flaps
- Smooth cannulation facilitated by tapered tip
- Adapts to bile duct due to gentle curved shape
- Pushability and softness delivered by polyethylene material
- Available in the widest range of sizes of any stent
- Simplify inventory management and enhance procedural efficiency with stent sets, which include pushing and guiding catheters




*Image courtesy of Dr Umesh M. Jalihal,  
M S Ramaiah Memorial Hospital, Bangalore, India*

## BILIARY STENTS

# Cotton-Leung<sup>®</sup>

## SOFT-FLEX<sup>®</sup> BILIARY STENT

The same subtle, natural curve as the Cotton-Leung Biliary Stent in a pliable material that is soft on tissue and strong on results.



### Cotton-Leung®

SOFT-FLEX® BILIARY STENT

- Fluoroscopic visibility enhanced by Sof-Flex, a soft pliable radiopaque material
- Minimizes migration with proximal and distal flaps
- Adapts to bile duct because of gentle curved shape
- Smooth cannulation facilitated by tapered tip
- Available in the widest range of sizes of any stent




## BILIARY STENTS

# Cotton-Huibregtse<sup>®</sup>

## BILIARY STENT

Shaping the future of biliary drainage, this plastic stent is distinguished by its anatomical design, which enhances positioning in the bile duct.



# Cotton-Huibregtse®

### BILIARY STENT

- Minimizes migration with proximal and distal flaps
- Smooth cannulation facilitated by tapered tip
- Enhanced positioning provided by anatomical shape
- Pushability and softness delivered by polyethylene material
- Simplify inventory management and enhance procedural efficiency with stent sets, which include pushing and guiding catheters




## BILIARY STENTS

# ST-2 Soehendra Tannenbaum®

### BILIARY STENT


The Teflon® material in our ST-2 Soehendra-Tannenbaum Stent is formulated to enhance patency with low surface friction to optimize flow.



# ST-2 Soehendra Tannenbaum®

### BILIARY STENT


- Minimizes stent migration with proximal and distal flap design
- Proper positioning in bile duct enabled by gentle contour
- Smoothness and pushability delivered by Teflon material
- Available in a wide range of sizes



# OASIS<sup>®</sup>

## ONE ACTION STENT INTRODUCTION SYSTEM - PRELOADED

Whether you choose a pre-loaded Cotton-Leung Stent or ST-2 Soehendra-Tannenbaum, Cook's OASIS One Action Stent Introduction System delivers fast, efficient stent placement that can streamline your plastic stenting procedures.



# OASIS®


### ONE ACTION STENT INTRODUCTION SYSTEM - PRELOADED

- Enhance efficiency and save time with single-step delivery
- Maximize efficiency with preloaded Cotton-Leung or ST-2 Soehendra-Tannenbaum stents
- Fluoroscopic visualization during placement is enhanced by radiopaque bands on introducer
- Stricture dilated prior to stent placement by tapered guiding catheter



*Image courtesy of Michael Kimmey, MD,  
Tacoma Digestive Disease Center, Tacoma, WA*

# Pancreatic Stents




## PANCREATIC STENTS

# Geenen<sup>®</sup>

### PANCREATIC STENT


With its contoured shape and offset flaps, the Geenen Pancreatic Stent is designed for optimal positioning with minimal migration.



### Geenen<sup>®</sup>

#### PANCREATIC STENT

- Smooth cannulation eased by tapered tip
- Proper positioning facilitated by contoured shape
- Drainage enhanced by spiral side holes
- Minimizes migration with offset flaps
- Pushability and softness delivered by the stent's polyethylene material
- Simplify inventory management and enhance procedure efficiency with stent sets, which include pushing catheter




*Image courtesy of Michael Gluck, MD,  
Virginia Mason Medical Center, Seattle, WA*

## PANCREATIC STENTS

# Geenen<sup>®</sup>

### SOF-FLEX<sup>®</sup> PANCREATIC STENT

The Geenen Sof-Flex Pancreatic Stent, made of a strong yet flexible material, is designed to gently take the shape of the duct instead of forcing the duct to take the shape of the stent. Available with and without internal ductal flaps to meet all of your clinical needs.



### Geenen<sup>®</sup>

#### SOF-FLEX<sup>®</sup> PANCREATIC STENT

- Fluoroscopic visibility enhanced by Sof-Flex, a soft pliable radiopaque material
- Smooth cannulation eased by tapered tip
- Proper positioning facilitated by contoured shape
- Drainage enhanced by spiral side holes
- Minimizes migration with offset flaps




*Image courtesy of Irving Waxman, MD,  
University of Chicago, Chicago, IL*

## PANCREATIC STENTS

# Zimmon<sup>®</sup>

### PANCREATIC STENTS

The pancreatic duct demands specialized care; Zimmon Pancreatic Stents are uniquely designed to meet these specific needs. The pigtail retention flap stent design aids in maintaining its position in the pancreatic duct.



### Zimmon®

#### PANCREATIC STENTS

- Minimize migration with duodenal pigtail (and optional ductal flap)
- Drainage enhanced by spiral side holes on select stents
- Smooth cannulation facilitated by tapered tip
- Accurate positioning enhanced by radiopaque markers on select stents
- Pushability and softness delivered by polyethylene material
- Meet specific procedural needs—with a variety of French sizes




## PANCREATIC STENTS

# Johlin™

### PANCREATIC WEDGE STENT

Designed to conform to the duct instead of forcing the duct to conform to the stent, Johlin Pancreatic Wedge Stents are made of flexible, yet strong Sof-Flex material. Available in a variety of sizes and can be smoothly delivered with the Wedge Stent Introducer.




## PANCREATIC STENTS

# Johlin™

PANCREATIC WEDGE STENT

- Eases entrance into duct and minimizes stent migration with wedge-shaped distal end
- Fluoroscopic visibility enhanced by Sof-Flex, a soft pliable radiopaque material
- Allows drainage of the lateral branches of the pancreatic duct through drainage holes
- Includes Wedge Stent Introducer, designed specifically for the Johlin Stent



*Image courtesy of Amrita Sethi, MD,  
Columbia University Medical Center, New York, NY*

# Other Stent Introducers




## INTRODUCERS

# Fusion<sup>®</sup>OASIS<sup>®</sup>

## ONE ACTION STENT INTRODUCTION SYSTEM


Streamline plastic stenting procedures with Fusion OASIS. Accurately reposition or remove plastic stents while maintaining a secure wire guide, reducing procedural steps and saving valuable time.



# Fusion® OASIS®

## ONE ACTION STENT INTRODUCTION SYSTEM

- Reposition, recapture before deployment or remove plastic stents entirely while maintaining ductal access thanks to Fusion capability
- Compatible with a wide array of our exceptional plastic stents.
- Reduce procedural steps by delivering multiple stents without recannulation
- Integrated guiding and pushing catheter provides optimal pushability when negotiating tortuous anatomy.




*Image courtesy of Kapil Gupta, MD,  
Hennepin County Medical Center, Minneapolis, MN*

## INTRODUCERS

# OASIS<sup>®</sup>

## ONE ACTION STENT INTRODUCTION SYSTEM

Maximize efficiency in plastic stent placement with this one step delivery system. Choose from a wide variety of our plastic stents to meet your clinical needs.



# OASIS®

### ONE ACTION STENT INTRODUCTION SYSTEM

- Enhance efficiency and save time with single-step delivery
- Fluoroscopic visualization during placement is enhanced by radiopaque bands on introducer
- Stricture dilated prior to stent placement by tapered guiding catheter



*Image courtesy of Dr Umesh M. Jalihal,  
M S Ramaiah Memorial Hospital, Bangalore, India*


# Introducers

Our variety of introducers gives you options for endoscopic stent placement.



# Introducers

- Pushing catheters available in a range of sizes and compatible with both straight and pigtail stents.
- Enhanced visualization delivered by radiopaque markers on our guiding catheters.
- Added convenience delivered by our stent introducer sets which include guiding catheter and pushing catheter components.



*Image courtesy of Dr. Rishi Pawa,  
Wake Forest Baptist Medical Center,  
Winston-Salem, NC*

# Order Information

|                          |      | 3 FR               | 4 FR               | 5 FR               | 6 FR               | 7 FR               | French Size                                                                                                             |
|--------------------------|------|--------------------|--------------------|--------------------|--------------------|--------------------|-------------------------------------------------------------------------------------------------------------------------|
| Length                   |      |                    |                    |                    |                    |                    |                                                                                                                         |
| Order Number             | 2 cm | GPSO-3-2<br>G27105 | GPSO-4-2<br>G27106 | GPSO-5-2<br>G22434 | GPSO-6-2<br>G22467 | GPSO-7-<br>G22379  | Color Connection<br>Purple accepts .035" wire guide<br>Black accepts .021" wire guide<br>Green accepts .018" wire guide |
|                          | 3 cm | GPSO-3-3<br>G27190 |                    | GPSO-5-3<br>G22107 | GPSO-6-3<br>G22410 | GPSO-7-3<br>G21371 |                                                                                                                         |
| Global Part Number (GPN) | 4 cm | GPSO-3-4<br>G24543 |                    | GPSO-5-4<br>G22110 |                    | GPSO-7-4<br>G21373 |                                                                                                                         |
|                          | 5    | GPSO-3-5           | GPSO-4-5           | GPSO-5-5<br>G22111 | GPSO-6-5<br>G22411 | GPSO-7-5<br>G21374 |                                                                                                                         |

BILIARY STENTS

Zimmon®  
BILIARY STENTS

|       | 5 FR               | 6 FR               | 7 FR               | 8 FR               | 10 FR               |
|-------|--------------------|--------------------|--------------------|--------------------|---------------------|
| 2 cm  |                    |                    |                    |                    | ZSO-10-2<br>G25399  |
| 3 cm  | ZSO-5-3<br>G22032  |                    | ZSO-7-3<br>G22159  |                    | ZSO-10-3<br>G21973  |
| 4 cm  | ZSO-5-4<br>G22034  | ZSO-6-4<br>G22105  | ZSO-7-4<br>G22161  | ZSO-8-4<br>G21369  | ZSO-10-4<br>G21974  |
| 5 cm  | ZSO-5-5<br>G22036  |                    | ZSO-7-5<br>G22163  | ZSO-8-5<br>G21370  | ZSO-10-5<br>G21975  |
| 6 cm  |                    |                    | ZSO-7-6<br>G22164  |                    | ZSO-10-6<br>G21976  |
| 7 cm  | ZSO-5-7<br>G22038  | ZSO-6-7<br>G22109  | ZSO-7-7<br>G22167  | ZSO-8-7<br>G21372  | ZSO-10-7<br>G21977  |
| 8 cm  |                    |                    |                    |                    | ZSO-10-8<br>G21980  |
| 9 cm  | ZSO-5-9<br>G22040  | ZSO-6-9<br>G22112  | ZSO-7-9<br>G22170  | ZSO-8-9<br>G21375  | ZSO-10-9<br>G21982  |
| 10 cm | ZSO-5-10<br>G21499 | ZSO-6-10<br>G21998 | ZSO-7-10<br>G21707 | ZSO-8-10<br>G21401 | ZSO-10-10<br>G21637 |
| 12 cm | ZSO-5-12<br>G21501 |                    | ZSO-7-12<br>G21708 | ZSO-8-12<br>G21404 | ZSO-10-12<br>G21640 |
| 15 cm | ZSO-5-15<br>G22443 |                    | ZSO-7-15<br>G21711 | ZSO-8-15<br>G21408 | ZSO-10-15<br>G21645 |
| 16 cm |                    |                    | ZSO-7-16<br>G23005 |                    |                     |
| 18 cm |                    |                    | ZSO-7-18<br>G23006 |                    |                     |

Components available separately.

BILIARY STENTS

Zimmon®  
BILIARY STENTS


|       | 5 FR                | 6 FR                | 7 FR                | 8 FR                | 10 FR                | 11.5 FR                |
|-------|---------------------|---------------------|---------------------|---------------------|----------------------|------------------------|
| 2 cm  |                     |                     | ZEBD-7-2<br>G22438  |                     | ZEBD-10-2<br>G22441  |                        |
| 3 cm  | ZEBD-5-3<br>G22436  |                     | ZEBD-7-3<br>G22439  |                     | ZEBD-10-3<br>G21780  |                        |
| 4 cm  | ZEBD-5-4<br>G22137  | ZEBD-6-4<br>G22191  | ZEBD-7-4<br>G21393  |                     | ZEBD-10-4<br>G21781  |                        |
| 5 cm  | ZEBD-5-5<br>G22138  |                     | ZEBD-7-5<br>G21394  | ZEBD-8-5<br>G21460  | ZEBD-10-5<br>G21783  | ZEBD-11.5-5<br>G21439  |
| 6 cm  |                     |                     | ZEBD-7-6<br>G21395  |                     | ZEBD-10-6<br>G21785  |                        |
| 7 cm  | ZEBD-5-7<br>G22140  | ZEBD-6-7<br>G22194  | ZEBD-7-7<br>G21396  | ZEBD-8-7<br>G21462  | ZEBD-10-7<br>G21787  | ZEBD-11.5-7<br>G21442  |
| 8 cm  |                     |                     | ZEBD-7-8<br>G21397  |                     | ZEBD-10-8<br>G21788  |                        |
| 9 cm  | ZEBD-5-9<br>G22141  |                     | ZEBD-7-9<br>G21398  | ZEBD-8-9<br>G22440  | ZEBD-10-9<br>G21789  | ZEBD-11.5-9<br>G21446  |
| 10 cm | ZEBD-5-10<br>G22122 | ZEBD-6-10<br>G21795 | ZEBD-7-10<br>G25794 |                     | ZEBD-10-10<br>G21862 |                        |
| 12 cm | ZEBD-5-12<br>G22437 | ZEBD-6-12<br>G21798 | ZEBD-7-12<br>G21526 | ZEBD-8-12<br>G22026 | ZEBD-10-12<br>G21864 | ZEBD-11.5-12<br>G22207 |
| 15 cm | ZEBD-5-15<br>G22128 |                     | ZEBD-7-15<br>G21529 | ZEBD-8-15<br>G22027 | ZEBD-10-15<br>G21865 | ZEBD-11.5-15<br>G22442 |
| 18 cm |                     |                     | ZEBD-7-18<br>G21530 |                     |                      |                        |

Set includes: Stent and 170 cm pushing catheter.


BILIARY STENTS

Zimmon®  
BILIARY STENTS

ZSO



ZEBD




Compass BDS™  
BILIARY PLASTIC STENT

|       |                     |
|-------|---------------------|
|       | 7 FR                |
| 5 cm  | CBSO-7-5<br>G55519  |
| 10 cm | CBSO-7-10<br>G55520 |
| 15 cm | CBSO-7-15<br>G55521 |

Components available separately.

CBSO




Solus<sup>®</sup>

DOUBLE PIGTAIL STENT WITH INTRODUCER

|       | 10 FR                  |
|-------|------------------------|
|       |                        |
| 1 cm  | ZSS-10-1-RB<br>G26829  |
| 3 cm  | ZSS-10-3-RB<br>G25670  |
| 4 cm  | ZSS-10-4-RB<br>G25671  |
| 5 cm  | ZSS-10-5-RB<br>G25672  |
| 7 cm  | ZSS-10-7-RB<br>G25673  |
| 9 cm  | ZSS-10-9-RB<br>G25674  |
| 10 cm | ZSS-10-10-RB<br>G25685 |
| 12 cm | ZSS-10-12-RB<br>G25675 |
| 15 cm | ZSS-10-15-RB<br>G25676 |

Set includes stent, pushing catheter, guiding catheter and pigtail straightener (if applicable). Recommend Tracer Metro<sup>®</sup>, sold separately.

ZSS



BILIARY STENTS

Cotton-Leung®  
BILIARY STENT

|       | 5 FR                | 7 FR                | 8.5 FR                | 10 FR                | 11.5 FR                |
|-------|---------------------|---------------------|-----------------------|----------------------|------------------------|
| 2 cm  |                     | CLSO-7-2<br>G22392  |                       | CLSO-10-2<br>G21463  |                        |
| 3 cm  | CLSO-5-3<br>G22383  | CLSO-7-3<br>G21679  | CLSO-8.5-3<br>G21617  | CLSO-10-3<br>G21464  | CLSO-11.5-3<br>G21639  |
| 4 cm  |                     | CLSO-7-4<br>G21680  |                       | CLSO-10-4<br>G21465  |                        |
| 5 cm  | CLSO-5-5<br>G21565  | CLSO-7-5<br>G21681  | CLSO-8.5-5<br>G21619  | CLSO-10-5<br>G22622  | CLSO-11.5-5<br>G21643  |
| 6 cm  |                     | CLSO-7-6<br>G27143  | CLSO-8.5-6<br>G27148  | CLSO-10-6<br>G21467  | CLSO-11.5-6<br>G27141  |
| 7 cm  | CLSO-5-7<br>G22384  | CLSO-7-7<br>G21683  | CLSO-8.5-7<br>G21620  | CLSO-10-7<br>G21469  | CLSO-11.5-7<br>G21646  |
| 8 cm  |                     | CLSO-7-8<br>G21684  | CLSO-8.5-8<br>G22343  | CLSO-10-8<br>G21472  | CLSO-11.5-8<br>G27104  |
| 9 cm  | CLSO-5-9<br>G22385  | CLSO-7-9<br>G21685  | CLSO-8.5-9<br>G21621  | CLSO-10-9<br>G21474  | CLSO-11.5-9<br>G21647  |
| 10 cm |                     | CLSO-7-10<br>G21979 | CLSO-8.5-10<br>G22009 | CLSO-10-10<br>G21807 | CLSO-11.5-10<br>G21761 |
| 11 cm |                     | CLSO-7-11<br>G21981 | CLSO-8.5-11<br>G27145 | CLSO-10-11<br>G21808 | CLSO-11.5-11<br>G27103 |
| 12 cm | CLSO-5-12<br>G22398 | CLSO-7-12<br>G21983 | CLSO-8.5-12<br>G22011 | CLSO-10-12<br>G21810 | CLSO-11.5-12<br>G27140 |
| 13 cm |                     | CLSO-7-13<br>G21984 | CLSO-8.5-13<br>G27146 | CLSO-10-13<br>G21811 | CLSO-11.5-13<br>G31508 |
| 14 cm |                     | CLSO-7-14<br>G21987 | CLSO-8.5-14<br>G27147 | CLSO-10-14<br>G21812 | CLSO-11.5-14<br>G21763 |
| 15 cm | CLSO-5-15<br>G22399 | CLSO-7-15<br>G21988 | CLSO-8.5-15<br>G22012 | CLSO-10-15<br>G21813 | CLSO-11.5-15<br>G21764 |
| 16 cm |                     | CLSO-7-16<br>G27142 | CLSO-8.5-16<br>G22709 | CLSO-10-16<br>G21815 | CLSO-11.5-16<br>G21766 |
| 17 cm |                     | CLSO-7-17<br>G22480 | CLSO-8.5-17<br>G22710 | CLSO-10-17<br>G21816 | CLSO-11.5-17<br>G22400 |

BILIARY STENTS

Cotton-Leung®  
BILIARY STENT

|       | 5 FR | 7 FR                | 8.5 FR                | 10 FR                | 11.5 FR                |
|-------|------|---------------------|-----------------------|----------------------|------------------------|
| 18 cm |      | CLSO-7-18<br>G21990 | CLSO-8.5-18<br>G22711 | CLSO-10-18<br>G21818 | CLSO-11.5-18<br>G21767 |
| 19 cm |      | CLSO-7-19<br>G21991 |                       | CLSO-10-19<br>G22479 |                        |
| 20 cm |      | CLSO-7-20<br>G21993 | CLSO-8.5-20<br>G22712 | CLSO-10-20<br>G21821 | CLSO-11.5-20<br>G22713 |
| 21 cm |      | CLSO-7-21<br>G24784 |                       | CLSO-10-21<br>G24785 |                        |

Components available separately.

|       | 5 FR                | 7 FR                | 8.5 FR                | 10 FR                | 11.5 FR                |
|-------|---------------------|---------------------|-----------------------|----------------------|------------------------|
| 3 cm  |                     | CLBS-7-3<br>G21947  | CLBS-8.5-3<br>G21487  | CLBS-10-3<br>G21841  | CLBS-11.5-3<br>G22051  |
| 4 cm  |                     | CLBS-7-4<br>G21949  |                       | CLBS-10-4<br>G21843  |                        |
| 5 cm  | CLBS-5-5<br>G22402  | CLBS-7-5<br>G21950  | CLBS-8.5-5<br>G21489  | CLBS-10-5<br>G21844  | CLBS-11.5-5<br>G22052  |
| 6 cm  |                     | CLBS-7-6<br>G21951  | CLBS-8.5-6<br>G21490  | CLBS-10-6<br>G27124  | CLBS-11.5-6<br>G22054  |
| 7 cm  | CLBS-5-7<br>G22278  | CLBS-7-7<br>G21952  | CLBS-8.5-7<br>G21492  | CLBS-10-7<br>G21847  | CLBS-11.5-7<br>G22056  |
| 8 cm  |                     | CLBS-7-8<br>G27135  | CLBS-8.5-8<br>G21495  | CLBS-10-8<br>G21848  | CLBS-11.5-8<br>G22057  |
| 9 cm  | CLBS-5-9<br>G22403  | CLBS-7-9<br>G21954  | CLBS-8.5-9<br>G21496  | CLBS-10-9<br>G21850  | CLBS-11.5-9<br>G22059  |
| 10 cm |                     | CLBS-7-10<br>G21563 | CLBS-8.5-10<br>G21598 | CLBS-10-10<br>G21668 | CLBS-11.5-10<br>G27125 |
| 11 cm |                     | CLBS-7-11<br>G27130 | CLBS-8.5-11<br>G27136 | CLBS-10-11<br>G27121 | CLBS-11.5-11<br>G22053 |
| 12 cm | CLBS-5-12<br>G22414 | CLBS-7-12<br>G21566 | CLBS-8.5-12<br>G21599 | CLBS-10-12<br>G21670 | CLBS-11.5-12<br>G22055 |
| 13 cm |                     | CLBS-7-13<br>G27131 | CLBS-8.5-13<br>G27137 | CLBS-10-13<br>G21672 | CLBS-11.5-13<br>G27126 |
| 14 cm |                     | CLBS-7-14<br>G27132 | CLBS-8.5-14<br>G27138 | CLBS-10-14<br>G27122 | CLBS-11.5-14<br>G27127 |


BILIARY STENTS

Cotton-Leung®  
BILIARY STENT


|       | 5 FR | 7 FR                | 8.5 FR                | 10 FR                | 11.5 FR                |
|-------|------|---------------------|-----------------------|----------------------|------------------------|
| 15 cm |      | CLBS-7-15<br>G21569 | CLBS-8.5-15<br>G21600 | CLBS-10-15<br>G21673 | CLBS-11.5-15<br>G22060 |
| 16 cm |      | CLBS-7-16<br>G27133 | CLBS-8.5-16<br>G24648 | CLBS-10-16<br>G27123 | CLBS-11.5-16<br>G27128 |
| 17 cm |      | CLBS-7-17<br>G27134 | CLBS-8.5-17<br>G24649 | CLBS-10-17<br>G21674 | CLBS-11.5-17<br>G27129 |
| 18 cm |      | CLBS-7-18<br>G21574 | CLBS-8.5-18<br>G27139 | CLBS-10-18<br>G21675 | CLBS-11.5-18<br>G22064 |

Set Includes: Stent, 170 cm pushing catheter and 320 cm guiding catheter. 5 and 7 FR stents do not include guiding catheter.

CLSO




CLBS



Cotton-Leung®  
SOF-FLEX® BILIARY STENT

|       | 7 FR                   | 10 FR                   |
|-------|------------------------|-------------------------|
| 5 cm  | CLSO-SF-7-5<br>G49658  | CLSO-SF-10-5<br>G49653  |
| 7 cm  | CLSO-SF-7-7<br>G49659  | CLSO-SF-10-7<br>G49654  |
| 9 cm  | CLSO-SF-7-9<br>G49660  | CLSO-SF-10-9<br>G49655  |
| 12 cm | CLSO-SF-7-12<br>G49661 | CLSO-SF-10-12<br>G49656 |
| 15 cm | CLSO-SF-7-15<br>G49662 | CLSO-SF-10-15<br>G49657 |

Components available separately.



BILIARY STENTS


Cotton-Huibregtse<sup>®</sup>

BILIARY STENT

|       | 7 FR                 | 8.5 FR                 | 10 FR                 | 11.5 FR                 |
|-------|----------------------|------------------------|-----------------------|-------------------------|
| 3 cm  |                      |                        | CHBSO-10-3<br>G21648  |                         |
| 4 cm  | CHBSO-7-4<br>G21719  |                        | CHBSO-10-4<br>G21650  |                         |
| 5 cm  | CHBSO-7-5<br>G21720  | CHBSO-8.5-5<br>G21699  | CHBSO-10-5<br>G21652  | CHBSO-11.5-5<br>G21797  |
| 6 cm  |                      |                        | CHBSO-10-6<br>G21654  |                         |
| 7 cm  | CHBSO-7-7<br>G21721  | CHBSO-8.5-7<br>G21700  | CHBSO-10-7<br>G21656  | CHBSO-11.5-7<br>G21799  |
| 8 cm  |                      |                        | CHBSO-10-8<br>G21657  | CHBSO-11.5-8<br>G21800  |
| 9 cm  | CHBSO-7-9<br>G21722  | CHBSO-8.5-9<br>G21702  | CHBSO-10-9<br>G21658  | CHBSO-11.5-9<br>G21801  |
| 10 cm | CHBSO-7-10<br>G22211 | CHBSO-8.5-10<br>G22195 | CHBSO-10-10<br>G21879 | CHBSO-11.5-10<br>G22197 |
| 11 cm |                      |                        | CHBSO-10-11<br>G21880 |                         |
| 12 cm | CHBSO-7-12<br>G22213 | CHBSO-8.5-12<br>G22196 | CHBSO-10-12<br>G21881 | CHBSO-11.5-12<br>G22199 |
| 15 cm | CHBSO-7-15<br>G22217 | CHBSO-8.5-15<br>G22198 | CHBSO-10-15<br>G21883 | CHBSO-11.5-15<br>G22200 |
| 18 cm | CHBSO-7-18<br>G27172 | CHBSO-8.5-18<br>G27173 | CHBSO-10-18<br>G22404 | CHBSO-11.5-18<br>G27171 |

Components available separately.


CHBSO



|       | 7 FR                | 8.5 FR                | 10 FR                | 11.5 FR                |
|-------|---------------------|-----------------------|----------------------|------------------------|
| 3 cm  | CHBS-7-3<br>G22486  |                       | CHBS-10-3<br>G21454  |                        |
| 5 cm  | CHBS-7-5<br>G21831  | CHBS-8.5-5<br>G21814  | CHBS-10-5<br>G21455  | CHBS-11.5-5<br>G21541  |
| 6 cm  | CHBS-7-6<br>G22485  |                       |                      |                        |
| 7 cm  | CHBS-7-7<br>G21832  | CHBS-8.5-7<br>G21817  | CHBS-10-7<br>G21457  | CHBS-11.5-7<br>G21542  |
| 9 cm  | CHBS-7-9<br>G21834  | CHBS-8.5-9<br>G21820  | CHBS-10-9<br>G21459  | CHBS-11.5-9<br>G21545  |
| 10 cm | CHBS-7-10<br>G22420 |                       | CHBS-10-10<br>G21996 | CHBS-11.5-10<br>G22118 |
| 11 cm |                     |                       | CHBS-10-11<br>G21997 |                        |
| 12 cm | CHBS-7-12<br>G21970 | CHBS-8.5-12<br>G21888 | CHBS-10-12<br>G21999 | CHBS-11.5-12<br>G22121 |
| 15 cm | CHBS-7-15<br>G21971 | CHBS-8.5-15<br>G21889 | CHBS-10-15<br>G22000 | CHBS-11.5-15<br>G22123 |

Set Includes: Stent, 170 cm pushing catheter and 320 cm guiding catheter. 7 FR stents do not include guiding catheter.

CHBS




ST-2 Soehendra Tannenbaum®

BILIARY STENT

|       | 8.5 FR                | 10 FR                | 11.5 FR                |
|-------|-----------------------|----------------------|------------------------|
| 5 cm  | TTSO-8.5-5<br>G22292  | TTSO-10-5<br>G21436  | TTSO-11.5-5<br>G22283  |
| 6 cm  | TTSO-8.5-6<br>G22293  | TTSO-10-6<br>G21437  | TTSO-11.5-6<br>G22284  |
| 7 cm  | TTSO-8.5-7<br>G22294  | TTSO-10-7<br>G21438  | TTSO-11.5-7<br>G22285  |
| 8 cm  | TTSO-8.5-8<br>G22295  | TTSO-10-8<br>G21440  | TTSO-11.5-8<br>G22286  |
| 9 cm  | TTSO-8.5-9<br>G22296  | TTSO-10-9<br>G21441  | TTSO-11.5-9<br>G22287  |
| 10 cm | TTSO-8.5-10<br>G22288 | TTSO-10-10<br>G22125 | TTSO-11.5-10<br>G22280 |
| 11 cm | TTSO-8.5-11<br>G22289 | TTSO-10-11<br>G22127 | TTSO-11.5-11<br>G27169 |
| 12 cm | TTSO-8.5-12<br>G22290 | TTSO-10-12<br>G22129 | TTSO-11.5-12<br>G22281 |
| 13 cm | TTSO-8.5-13<br>G27198 | TTSO-10-13<br>G27167 | TTSO-11.5-13<br>G27196 |
| 14 cm | TTSO-8.5-14<br>G27170 | TTSO-10-14<br>G27168 | TTSO-11.5-14<br>G27197 |
| 15 cm | TTSO-8.5-15<br>G22291 | TTSO-10-15<br>G22130 | TTSO-11.5-15<br>G22282 |

Components available separately. Recommended for use with OASIS Introduction System.

TTSO




BILIARY STENTS

OASIS®  
ONE ACTION STENT INTRODUCTION SYSTEM

Preloaded with Cotton-Leung  
Biliary Stents

|       | 7 FR                | 8.5 FR                | 10 FR                | 11.5 FR                |
|-------|---------------------|-----------------------|----------------------|------------------------|
| 5 cm  | OACL-7-5<br>G25656  | OACL-8.5-5<br>G25374  | OACL-10-5<br>G25379  | OACL-11.5-5<br>G25384  |
| 7 cm  | OACL-7-7<br>G25657  | OACL-8.5-7<br>G25375  | OACL-10-7<br>G25380  | OACL-11.5-7<br>G25385  |
| 9 cm  | OACL-7-9<br>G25658  | OACL-8.5-9<br>G25376  | OACL-10-9<br>G25381  | OACL-11.5-9<br>G25386  |
| 12 cm | OACL-7-12<br>G25659 | OACL-8.5-12<br>G25377 | OACL-10-12<br>G25382 | OACL-11.5-12<br>G25387 |
| 15 cm | OACL-7-15<br>G25660 | OACL-8.5-15<br>G25378 | OACL-10-15<br>G25383 | OACL-11.5-15<br>G25388 |


Set includes: introducer with four radiopaque bands 5 cm apart, preloaded Cotton-Leung stent, 170 cm pushing catheter and 205 cm guiding catheter. Recommended Tracer Metro® wire guide, sold separately.



Preloaded with ST-2 Soehendra-Tannenbaum  
Biliary Stents

|       | 8.5 FR                | 10 FR                | 11.5 FR                |
|-------|-----------------------|----------------------|------------------------|
| 5 cm  | OATS-8.5-5<br>G26846  | OATS-10-5<br>G26858  | OATS-11.5-5<br>G26836  |
| 7 cm  | OATS-8.5-7<br>G26848  | OATS-10-7<br>G26859  | OATS-11.5-7<br>G26837  |
| 9 cm  | OATS-8.5-9<br>G26849  | OATS-10-9<br>G26860  | OATS-11.5-9<br>G26838  |
| 12 cm | OATS-8.5-12<br>G26853 | OATS-10-12<br>G26863 | OATS-11.5-12<br>G26842 |
| 15 cm | OATS-8.5-15<br>G26854 | OATS-10-15<br>G26864 | OATS-11.5-15<br>G26843 |

Set includes: introducer with four radiopaque bands 5 cm apart, preloaded ST-2 Soehendra-Tannenbaum stent, 170 cm pushing catheter and 205 cm guiding catheter. Recommended Tracer Metro® wire guide, sold separately.



PANCREATIC STENTS

Geenen®  
PANCREATIC STENT

|       | 3 FR                | 4 FR               | 5 FR                | 6 FR                | 7 FR                | 8.5 FR                | 10 FR                | 11.5 FR                |
|-------|---------------------|--------------------|---------------------|---------------------|---------------------|-----------------------|----------------------|------------------------|
| 2 cm  | GPSO-3-2<br>G27105  | GPSO-4-2<br>G27106 | GPSO-5-2<br>G22434  | GPSO-6-2<br>G22467  | GPSO-7-2<br>G22379  |                       | GPSO-10-2<br>G22465  |                        |
| 3 cm  | GPSO-3-3<br>G27190  |                    | GPSO-5-3<br>G22107  | GPSO-6-3<br>G22410  | GPSO-7-3<br>G21371  | GPSO-8.5-3<br>G22314  | GPSO-10-3<br>G22319  | GPSO-11.5-3<br>G22464  |
| 4 cm  | GPSO-3-4<br>G24543  |                    | GPSO-5-4<br>G22110  |                     | GPSO-7-4<br>G21373  |                       |                      |                        |
| 5 cm  | GPSO-3-5<br>G24545  | GPSO-4-5<br>G24525 | GPSO-5-5<br>G22111  | GPSO-6-5<br>G22411  | GPSO-7-5<br>G21374  | GPSO-8.5-5<br>G22315  | GPSO-10-5<br>G22275  | GPSO-11.5-5<br>G22389  |
| 6 cm  | GPSO-3-6<br>G27191  |                    | GPSO-5-6<br>G27162  |                     | GPSO-7-6<br>G21376  |                       |                      |                        |
| 7 cm  | GPSO-3-7<br>G24523  | GPSO-4-7<br>G24529 | GPSO-5-7<br>G22113  | GPSO-6-7<br>G22412  | GPSO-7-7<br>G21377  | GPSO-8.5-7<br>G22316  | GPSO-10-7<br>G22320  | GPSO-11.5-7<br>G22463  |
| 8 cm  | GPSO-3-8<br>G27192  |                    | GPSO-5-8<br>G27163  |                     | GPSO-7-8<br>G27166  |                       |                      |                        |
| 9 cm  | GPSO-3-9<br>G27193  | GPSO-4-9<br>G24801 | GPSO-5-9<br>G22114  | GPSO-6-9<br>G22413  | GPSO-7-9<br>G21379  | GPSO-8.5-9<br>G22317  | GPSO-10-9<br>G22322  | GPSO-11.5-9<br>G22462  |
| 10 cm | GPSO-3-10<br>G24542 |                    | GPSO-5-10<br>G22097 |                     | GPSO-7-10<br>G21506 |                       | GPSO-10-10<br>G22321 |                        |
| 11 cm | GPSO-3-11<br>G27185 |                    | GPSO-5-11<br>G27159 |                     | GPSO-7-11<br>G22466 |                       |                      |                        |
| 12 cm | GPSO-3-12<br>G27186 |                    | GPSO-5-12<br>G22099 | GPSO-6-12<br>G22401 | GPSO-7-12<br>G21507 | GPSO-8.5-12<br>G22318 | GPSO-10-12<br>G22323 |                        |
| 13 cm | GPSO-3-13<br>G27187 |                    | GPSO-5-13<br>G27160 |                     | GPSO-7-13<br>G27164 |                       |                      |                        |
| 14 cm | GPSO-3-14<br>G27188 |                    | GPSO-5-14<br>G27161 |                     | GPSO-7-14<br>G27165 |                       |                      |                        |
| 15 cm | GPSO-3-15<br>G27189 |                    | GPSO-5-15<br>G22102 | GPSO-6-15<br>G24546 | GPSO-7-15<br>G21509 | GPSO-8.5-15<br>G22986 | GPSO-10-15<br>G22987 | GPSO-11.5-15<br>G24541 |
| 18 cm |                     |                    | GPSO-5-18<br>G24802 |                     | GPSO-7-18<br>G22984 |                       |                      |                        |
| 20 cm |                     |                    |                     |                     | GPSO-7-20<br>G22985 |                       |                      |                        |

Components available separately. Recommend .018 Roadrunner® wire guide for 3 FR stents. Recommend Tracer Metro® Direct for all others.

PANCREATIC STENTS

Geenen®  
PANCREATIC STENT

|       | 3 FR                | 5 FR                | 7 FR                | 8.5 FR                | 10 FR                | 11.5 FR                |
|-------|---------------------|---------------------|---------------------|-----------------------|----------------------|------------------------|
| 2 cm  |                     | GEPD-5-2<br>G22478  |                     |                       |                      |                        |
| 3 cm  | GEPD-3-3<br>G24552  | GEPD-5-3<br>G22206  | GEPD-7-3<br>G21475  | GEPD-8.5-3<br>G22475  | GEPD-10-3<br>G22276  |                        |
| 4 cm  | GEPD-3-4<br>G27180  | GEPD-5-4<br>G22208  | GEPD-7-4<br>G21476  |                       |                      |                        |
| 5 cm  | GEPD-3-5<br>G27181  | GEPD-5-5<br>G22209  | GEPD-7-5<br>G21478  | GEPD-8.5-5<br>G22277  | GEPD-10-5<br>G21851  |                        |
| 6 cm  | GEPD-3-6<br>G27149  | GEPD-5-6<br>G22477  | GEPD-7-6<br>G27157  |                       |                      | GEPD-11.5-6<br>G22470  |
| 7 cm  | GEPD-3-7<br>G27182  | GEPD-5-7<br>G22210  | GEPD-7-7<br>G21479  | GEPD-8.5-7<br>G22474  | GEPD-10-7<br>G22472  |                        |
| 8 cm  | GEPD-3-8<br>G27183  | GEPD-5-8<br>G22476  | GEPD-7-8<br>G27158  |                       |                      | GEPD-11.5-8<br>G22469  |
| 9 cm  | GEPD-3-9<br>G27184  | GEPD-5-9<br>G22212  | GEPD-7-9<br>G21480  | GEPD-8.5-9<br>G22473  | GEPD-10-9<br>G22471  |                        |
| 10 cm | GEPD-3-10<br>G27174 | GEPD-5-10<br>G27150 | GEPD-7-10<br>G21572 |                       |                      | GEPD-11.5-10<br>G22468 |
| 11 cm | GEPD-3-11<br>G27175 | GEPD-5-11<br>G27151 | GEPD-7-11<br>G27154 |                       |                      |                        |
| 12 cm | GEPD-3-12<br>G27176 | GEPD-5-12<br>G22189 | GEPD-7-12<br>G21575 | GEPD-8.5-12<br>G24540 | GEPD-10-12<br>G22487 | GEPD-11.5-12<br>G22386 |
| 13 cm | GEPD-3-13<br>G27177 | GEPD-5-13<br>G27152 | GEPD-7-13<br>G27155 |                       |                      |                        |
| 14 cm | GEPD-3-14<br>G27178 | GEPD-5-14<br>G27153 | GEPD-7-14<br>G27156 |                       |                      |                        |
| 15 cm | GEPD-3-15<br>G27179 | GEPD-5-15<br>G22190 | GEPD-7-15<br>G21578 |                       | GEPD-10-15<br>G22488 |                        |

Set Includes: Stent, 170 cm pushing catheter and 320 cm guiding catheter. Recommend .018 Roadrunner® wire guide for 3 FR stents. Recommend Tracer Metro® Direct for all others.

PANCREATIC STENTS

Geenen®  
PANCREATIC STENT

|       | 5 FR                 | 7 FR                 |
|-------|----------------------|----------------------|
| 2 cm  | GPSOS-5-2<br>G22460  | GPSOS-7-2<br>G22459  |
| 3 cm  | GPSOS-5-3<br>G21513  | GPSOS-7-3<br>G21592  |
| 4 cm  | GPSOS-5-4<br>G21514  |                      |
| 5 cm  | GPSOS-5-5<br>G21515  | GPSOS-7-5<br>G21594  |
| 7 cm  | GPSOS-5-7<br>G21517  | GPSOS-7-7<br>G21595  |
| 9 cm  | GPSOS-5-9<br>G21519  | GPSOS-7-9<br>G21597  |
| 10 cm |                      | GPSOS-7-10<br>G21856 |
| 12 cm | GPSOS-5-12<br>G21633 | GPSOS-7-12<br>G21858 |
| 15 cm | GPSOS-5-15<br>G21636 | GPSOS-7-15<br>G25073 |

|       | 5 FR                | 7 FR                |
|-------|---------------------|---------------------|
| 2 cm  | GPDS-5-2<br>G22461  |                     |
| 3 cm  | GPDS-5-3<br>G21802  | GPDS-7-3<br>G21898  |
| 5 cm  | GPDS-5-5<br>G21803  | GPDS-7-5<br>G21899  |
| 7 cm  | GPDS-5-7<br>G21804  | GPDS-7-7<br>G21901  |
| 9 cm  | GPDS-5-9<br>G21806  | GPDS-7-9<br>G21902  |
| 12 cm | GPDS-5-12<br>G21723 | GPDS-7-12<br>G21939 |
| 15 cm | GPDS-5-15<br>G21724 |                     |


Set Includes: Stent, 170 cm pushing catheter  
and 320 cm guiding catheter. Recommend Tracer  
Metro® Direct wire guide.

Components available separately.


PANCREATIC STENTS

Geenen®  
PANCREATIC STENT


GPSO




GPSOS



GEPD



GPDS



PANCREATIC STENTS

Geenen®  
SOF-FLEX® PANCREATIC STENT


|       | 5 FR                   |
|-------|------------------------|
| 3 cm  | GPSO-SF-5-3<br>G49663  |
| 5 cm  | GPSO-SF-5-5<br>G49664  |
| 7 cm  | GPSO-SF-5-7<br>G49665  |
| 9 cm  | GPSO-SF-5-9<br>G49666  |
| 12 cm | GPSO-SF-5-12<br>G49667 |

Components available separately. Recommend Tracer Metro® Direct, sold separately.


|       | 5 FR                    |
|-------|-------------------------|
| 3 cm  | GPSOS-SF-5-3<br>G49668  |
| 5 cm  | GPSOS-SF-5-5<br>G49669  |
| 7 cm  | GPSOS-SF-5-7<br>G49670  |
| 9 cm  | GPSOS-SF-5-9<br>G49671  |
| 12 cm | GPSOS-SF-5-12<br>G49672 |

Components available separately. Recommend Tracer Metro® Direct, sold separately.

GPSO-SF



GPSOS-SF



PANCREATIC STENTS



|        | 4 FR                 | 5 FR                  | 6 FR                  | 7 FR                  | 8.5 FR                 | 10 FR                 |
|--------|----------------------|-----------------------|-----------------------|-----------------------|------------------------|-----------------------|
| 2 cm   | SPSOF-4-2<br>G24571  | SPSOF-5-2<br>G22357   | SPSOF-6-2<br>G22364   | SPSOF-7-2<br>G22350   |                        |                       |
| 2.5 cm |                      | SPSOF-5-2.5<br>G22376 | SPSOF-6-2.5<br>G22377 | SPSOF-7-2.5<br>G22378 |                        |                       |
| 3 cm   |                      | SPSOF-5-3<br>G22979   |                       | SPSOF-7-3<br>G22454   |                        | SPSOF-10-3<br>G22975  |
| 4 cm   | SPSOF-4-4<br>G24575  | SPSOF-5-4<br>G22358   | SPSOF-6-4<br>G22365   | SPSOF-7-4<br>G22351   |                        | SPSOF-10-4<br>G22976  |
| 5 cm   |                      | SPSOF-5-5<br>G22980   |                       | SPSOF-7-5<br>G22455   |                        | SPSOF-10-5<br>G24589  |
| 6 cm   | SPSOF-4-6<br>G24579  | SPSOF-5-6<br>G22359   | SPSOF-6-6<br>G22366   | SPSOF-7-6<br>G22352   |                        | SPSOF-10-6<br>G22977  |
| 7 cm   |                      | SPSOF-5-7<br>G22360   |                       | SPSOF-7-7<br>G22353   | SPSOF-8.5-7<br>G24590  | SPSOF-10-7<br>G24591  |
| 8 cm   | SPSOF-4-8<br>G24582  | SPSOF-5-8<br>G22361   | SPSOF-6-8<br>G22367   | SPSOF-7-8<br>G22354   |                        | SPSOF-10-8<br>G24707  |
| 9 cm   |                      | SPSOF-5-9<br>G22981   |                       | SPSOF-7-9<br>G22456   |                        |                       |
| 10 cm  | SPSOF-4-10<br>G24569 | SPSOF-5-10<br>G22362  | SPSOF-6-10<br>G22368  | SPSOF-7-10<br>G22355  | SPSOF-8.5-10<br>G24593 | SPSOF-10-10<br>G22982 |
| 11 cm  |                      | SPSOF-5-11<br>G27194  |                       | SPSOF-7-11<br>G27195  |                        |                       |
| 12 cm  | SPSOF-4-12<br>G50281 | SPSOF-5-12<br>G22363  | SPSOF-6-12<br>G22369  | SPSOF-7-12<br>G22356  | SPSOF-8.5-12<br>G24594 | SPSOF-10-12<br>G24592 |
| 15 cm  | SPSOF-4-15<br>G50282 | SPSOF-5-15<br>G22978  | SPSOF-6-15<br>G24584  | SPSOF-7-15<br>G22983  |                        | SPSOF-10-15<br>G24585 |
| 18 cm  | SPSOF-4-18<br>G26839 |                       |                       | SPSOF-7-18<br>G24586  |                        |                       |

Components available separately.

PANCREATIC STENTS

Zimmon®  
PANCREATIC STENTS

|       | 3 FR                   | 5 FR                | 7 FR                |
|-------|------------------------|---------------------|---------------------|
| 3 cm  |                        | SPSOS-5-3<br>G24570 | SPSOS-7-3<br>G24587 |
| 4 cm  | SPSOS-3-4-N<br>G24566  |                     |                     |
| 5 cm  |                        | SPSOS-5-5<br>G24574 | SPSOS-7-5<br>G24588 |
| 6 cm  | SPSOS-3-6-N<br>G24567  |                     |                     |
| 7 cm  |                        | SPSOS-5-7<br>G24578 |                     |
| 8 cm  | SPSOS-3-8-N<br>G24568  |                     |                     |
| 10 cm | SPSOS-3-10-N<br>G24564 |                     |                     |
| 12 cm | SPSOS-3-12-N<br>G24565 |                     |                     |
| 18 cm | SPSOS-3-18-N<br>G26847 |                     |                     |

Components available separately.

|       | 5 FR                 | 7 FR                 |
|-------|----------------------|----------------------|
| 2 cm  | ZPSOF-5-2<br>G22071  | ZPSOF-7-2<br>G22184  |
| 3 cm  | ZPSOF-5-3<br>G22375  |                      |
| 4 cm  | ZPSOF-5-4<br>G22073  | ZPSOF-7-4<br>G22186  |
| 5 cm  | ZPSOF-5-5<br>G22448  |                      |
| 6 cm  | ZPSOF-5-6<br>G22075  | ZPSOF-7-6<br>G22187  |
| 7 cm  | ZPSOF-5-7<br>G22076  | ZPSOF-7-7<br>G22451  |
| 8 cm  | ZPSOF-5-8<br>G22077  | ZPSOF-7-8<br>G22188  |
| 9 cm  | ZPSOF-5-9<br>G22449  | ZPSOF-7-9<br>G22452  |
| 10 cm | ZPSOF-5-10<br>G21827 | ZPSOF-7-10<br>G22065 |
| 12 cm | ZPSOF-5-12<br>G21828 | ZPSOF-7-12<br>G22068 |
| 15 cm | ZPSOF-5-15<br>G22450 | ZPSOF-7-15<br>G22435 |

Components available separately.

|       | 5 FR                 | 7 FR                    |
|-------|----------------------|-------------------------|
| 3 cm  | ZPSOS-5-3<br>G22453  |                         |
| 4 cm  |                      | ZPSOS-7-4<br>G21953     |
| 5 cm  | ZPSOS-5-5<br>G22457  |                         |
| 6 cm  |                      | ZPSOS-7-6<br>G21955     |
| 7 cm  | ZPSOS-5-7<br>G21849  |                         |
| 8 cm  |                      | ZPSOS-7-8<br>G21956     |
| 9 cm  | ZPSOS-5-9<br>G22458  |                         |
| 10 cm |                      | ZPSOS-7-10<br>G21638    |
| 12 cm | ZPSOS-5-12<br>G21419 | ZPSOS-7-12<br>G21642    |
| 25 cm |                      | ZPSOS-7-25-NP<br>G24610 |

Components available separately.

PANCREATIC STENTS

Zimmon®  
PANCREATIC STENTS

|       | 5 FR                 | 7 FR                 |
|-------|----------------------|----------------------|
| 2 cm  | ZEPDF-5-2<br>G22162  | ZEPDF-7-2<br>G21425  |
| 4 cm  | ZEPDF-5-4<br>G22165  | ZEPDF-7-4<br>G21426  |
| 5 cm  | ZEPDF-5-5<br>G22444  | ZEPDF-7-5<br>G22446  |
| 6 cm  | ZEPDF-5-6<br>G22169  | ZEPDF-7-6<br>G21427  |
| 7 cm  | ZEPDF-5-7<br>G22171  | ZEPDF-7-7<br>G22447  |
| 8 cm  | ZEPDF-5-8<br>G22173  | ZEPDF-7-8<br>G21428  |
| 9 cm  | ZEPDF-5-9<br>G22445  | ZEPDF-7-9<br>G21429  |
| 10 cm | ZEPDF-5-10<br>G22115 | ZEPDF-7-10<br>G21522 |
| 12 cm | ZEPDF-5-12<br>G22116 | ZEPDF-7-12<br>G21524 |
| 15 cm |                      | ZEPDF-7-15<br>G24563 |

Set Includes: Stent and 170 cm pushing catheter.


|       | 5 FR                 | 7 FR                 |
|-------|----------------------|----------------------|
| 2 cm  | ZEPDS-5-2<br>G21920  |                      |
| 4 cm  | ZEPDS-5-4<br>G21923  | ZEPDS-7-4<br>G22041  |
| 6 cm  | ZEPDS-5-6<br>G21924  | ZEPDS-7-6<br>G22042  |
| 8 cm  | ZEPDS-5-8<br>G21926  | ZEPDS-7-8<br>G22045  |
| 10 cm | ZEPDS-5-10<br>G21689 | ZEPDS-7-10<br>G21893 |
| 12 cm | ZEPDS-5-12<br>G21691 | ZEPDS-7-12<br>G21895 |

Set Includes: Stent and 170 cm pushing catheter.


PANCREATIC STENTS

Zimmon®  
PANCREATIC STENTS


SPSOF




ZPSOF




ZEPDF




SPSOS



ZPSOS



ZEPDS





PANCREATIC STENTS

Johlin™  
PANCREATIC WEDGE STENT

|       | 8.5 FR                | 10 FR                |
|-------|-----------------------|----------------------|
| 8 cm  | JPWS-8.5-8<br>G25694  | JPWS-10-8<br>G25713  |
| 10 cm | JPWS-8.5-10<br>G25695 | JPWS-10-10<br>G25714 |
| 12 cm | JPWS-8.5-12<br>G25693 | JPWS-10-12<br>G25715 |
| 14 cm | JPWS-8.5-14<br>G25691 | JPWS-10-14<br>G25716 |
| 16 cm | JPWS-8.5-16<br>G25690 | JPWS-10-16<br>G25717 |
| 18 cm | JPWS-8.5-18<br>G26830 | JPWS-10-18<br>G26825 |
| 20 cm | JPWS-8.5-20<br>G26831 | JPWS-10-20<br>G26827 |
| 22 cm | JPWS-8.5-22<br>G26832 | JPWS-10-22<br>G26828 |

JPWS



INTRODUCERS


Fusion<sup>®</sup>OASIS<sup>®</sup>  
ONE ACTION STENT INTRODUCTION SYSTEM

Fusion OASIS

|  | 8.5 FR              | 10 FR              |
|--|---------------------|--------------------|
|  |                     |                    |
|  | FS-OA-8.5<br>G31526 | FS-OA-10<br>G31527 |

170 cm pushing catheter and 195 cm guiding catheter. 4.2 mm minimum accessory channel.  
Not for use with pigtail stents.

FS-OA




Fusion Pushing Catheter

|  | 5 FR              | 7 FR              |
|--|-------------------|-------------------|
|  |                   |                   |
|  | FS-PC-5<br>G31550 | FS-PC-7<br>G31551 |

170 cm pushing catheter. 3.2 mm minimum accessory channel.

FS-PC



Pushing Catheter


|        | 3 FR           | 5 FR           | 7 FR           | 8.5 FR           | 10 FR           | 11.5 FR           |
|--------|----------------|----------------|----------------|------------------|-----------------|-------------------|
|        |                |                |                |                  |                 |                   |
| 170 cm | PC-3<br>G25092 | PC-5<br>G21771 | PC-7<br>G21774 | PC-8.5<br>G21623 | PC-10<br>G21443 | PC-11.5<br>G22091 |

PC



Guiding Catheter

|        | 5 FR           | 6 FR           |
|--------|----------------|----------------|
| 320 cm | GC-5<br>G22004 | GC-6<br>G22005 |



Huibregtse Guiding Catheter


|        | 6 FR            |
|--------|-----------------|
| 220 cm | HGC-6<br>G21448 |



Stent Introducer Sets

|        | 8.5 FR            | 10 FR            | 11.5 FR            |
|--------|-------------------|------------------|--------------------|
| 320 cm | SIS-8.5<br>G21532 | SIS-10<br>G22684 | SIS-11.5<br>G21749 |
| 170 cm |                   |                  |                    |

SIS





## Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)

Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)

Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)

Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)

Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)

Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)

France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)

Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)

Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)

Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)

Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)

Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)

Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)

Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)

Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)

Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)

Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)

Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)

United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

## Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)