

Streamline all your stone management procedures.



Fusion[®] Quattro[®]

EXTRACTION BALLOON

DUAL PLATFORM ERCP

ACCESS

STONE MANAGEMENT

STRICTURE MANAGEMENT



www.cookmedical.com



Fusion® Quattro®

EXTRACTION BALLOON

Used for endoscopic removal of stones in the biliary system and for contrast injection.

- 4 inflation diameters in 1 device to meet your specific clinical needs
- Ease of injection with contrast exiting the distal tip for above injection balloons
- Radiopaque bands enhance fluoroscopic visualization of the extraction balloon
- Reference marks on syringe indicate each inflation diameter

Order Number	Reference Part Number	Balloon Diameter mm	Catheter Size Fr	Catheter Length cm	Injection Port	Wire Guide Diameter inch
G48225	FS-QEB-A	8.5/10/12/15	6.6	200	Above Balloon	.035
G48226	FS-QEB-B	8.5/10/12/15	6.6	200	Below Balloon	.035
G31921	FS-QEB-XL-A	12/15/18/20	6.6	200	Above Balloon	.035
G31922	FS-QEB-XL-B	12/15/18/20	6.6	200	Below Balloon	.035

Minimum accessory channel 3.2 mm. **Caution:** This product contains natural rubber latex, which may cause an allergic reaction.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 97251996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com

