


ORDERING INFORMATION



Used to endoscopically dilate strictures of the esophagus.

Supplied sterile and is disposable - intended for single use only.

- Rapid balloon deflation* is achieved with the Express Evacuation Catheter and patent-pending Rapid Deflation Sleeve
- Kink-resistant Nitinol catheter provides maneuverability
- Flexible atraumatic tip for smooth access through especially tight strictures
- Excellent through-the-balloon visualization permits view of dilation progress
- Glow-in-the-dark size and inflation information on packaging and catheter tags for easy identification
- Pre-lubricated balloon facilitates removal from the accessory channel

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Inflated Balloon Diameter Fr	Balloon Length cm	ATM Pressure	PSI Pressure
G31928	HBD-8-9-10	8-9-10	24-27-30	8	4-6-8	60-90-120
G31925	HBD-10-11-12	10-11-12	30-33-36	8	2-4-6	30-60-90
G31926	HBD-12-13.5-15	12-13.5-15	36-40.5-45	8	2-4-6	30-60-90
G48732	HBD-15-16.5-18	15-16.5-18	45-49.5-54	8	2-4-6	30-60-90
G31927	HBD-18-19-20	18-19-20	54-57-60	8	2-4-6	30-60-90

Note: Catheter size is 7 Fr. Catheter length 180 cm. Minimum accessory channel 2.8 mm.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

*Based on using a 60cc syringe. Data on file at Cook Medical.

Customer Service Centers

AORTIC INTERVENTION
CRITICAL CARE
ENDOSCOPY
INTERVENTIONAL RADIOLOGY
LEAD MANAGEMENT
PERIPHERAL INTERVENTION
SURGERY
UROLOGY
WOMEN'S HEALTH

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com