

Deliver a fine staining-solution mist.



Glo-Tip®
SPRAY CATHETER

When spraying the gastrointestinal mucosa to identify suspicious lesions, a fine mist is essential to optimally stain the tissue. The Glo-Tip Spray Catheter delivers an exceptionally fine mist, allowing you to perform effective procedures. The 7 FR Teflon® sheath provides low coefficient of friction for even flow of the staining solution and for ease of insertion.



www.cookmedical.com



Glo-Tip®

SPRAY CATHETER

Used in conjunction with an endoscope to spray the gastrointestinal mucosa.

Supplied sterile and is disposable – intended for single use only.

- 7 Fr Teflon® sheath provides low coefficient to friction for ease of insertion, allows for even flow of staining solution
- Radiopaque tip enhances fluoroscopic visualization, if needed
- Smooth metal spray tip is atraumatic to mucosa
- Disposability eliminates cleaning and reprocessing costs, saves time

Order Number	Reference Part Number	French Size Fr	Catheter Length cm	Stylet Length cm	Tip
G24892	GT-7-SPRAY	7	240	100	Spray

Minimum accessory channel 2.8 mm.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI – www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 972519996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI – www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: cau.custserv@cookmedical.com

