


Tailor your treatment to the shape of the pancreatic duct.



Johlin™
PANCREATIC WEDGE STENT

Rather than forcing the pancreatic duct to conform to the shape of the stent, the Johlin conforms to the delicate shape of the duct. That's because we designed it with a wedge-shaped distal end to facilitate entrance into the duct and minimize stent migration. Made of flexible yet strong Sof-Flex® material, these stents are available in a wide variety of sizes to meet varying procedural needs and are smoothly delivered with the Wedge Stent Introducer.



www.cookmedical.com



Johlin™

PANCREATIC WEDGE STENT

Used to drain obstructed pancreatic ducts.

Supplied sterile and is disposable – intended for single use only.

- Made of soft, pliable Sof-Flex, a radiopaque material designed for strength and flexibility
- Engineered to conform to shape of pancreatic duct
- Ample drainage holes that potentially allow lateral branches of pancreatic duct to empty efficiently
- Packaged with its own introducer

Order Number	Reference Part Number	Stent French Size Fr	Stent Length cm	Introducer	Color Connection
G25694	JPWS-8.5-8	8.5	8	Preloaded	Purple
G25695	JPWS-8.5-10	8.5	10	Preloaded	Purple
G25693	JPWS-8.5-12	8.5	12	Preloaded	Purple
G25691	JPWS-8.5-14	8.5	14	Preloaded	Purple
G25690	JPWS-8.5-16	8.5	16	Preloaded	Purple
G26830	JPWS-8.5-18	8.5	18	Preloaded	Purple
G26831	JPWS-8.5-20	8.5	20	Preloaded	Purple
G26832	JPWS-8.5-22	8.5	22	Preloaded	Purple
G25713	JPWS-10-8	10	8	Preloaded	Purple
G25714	JPWS-10-10	10	10	Preloaded	Purple
G25715	JPWS-10-12	10	12	Preloaded	Purple
G25716	JPWS-10-14	10	14	Preloaded	Purple
G25717	JPWS-10-16	10	16	Preloaded	Purple
G26825	JPWS-10-18	10	18	Preloaded	Purple
G26827	JPWS-10-20	10	20	Preloaded	Purple
G26828	JPWS-10-22	10	22	Preloaded	Purple

Length measured from distal end to proximal end. Recommend Tracer Metro® .035 inch wire guide, sold separately.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 97251996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: cau.custserv@cookmedical.com

