Multiple solutions for scope cleaning.



Disposable Endoscopic Cleaning Brushes

To meet your varying endoscope-cleaning needs, Cook offers a wide array of Disposable Endoscopic Cleaning Brushes designed for cleaning accessory channels of upper and lower scopes. These brushes are available in double-ended configurations that allow double brushing with one pass and all feature covered brush tips designed to protect scope channels.





Disposable Endoscopic Cleaning Brushes

Used for cleaning the accessory channels of endoscopes.

Supplied non-sterile and is disposable - intended for single use only.

Endoscope valves, valve-housings and connectors can be cleaned using the Disposable Valve Cleaning Brush.

- Multiple configurations available to suit cleaning preferences
- Rigid, yet flexible channel brush shaft to allow for ease of passage
- Covered brush tip to protect endoscope channel
- Ergonomically designed double-ended valve cleaning brush to maximise cleaning

Order Number	Reference Part Number	Channel Brush	Brush Size mm	Valve Brush Size mm	Brush Length cm	Quantity Per Box
G30651	DCB-S-50	Single-Ended	6.5-5	n/a	240	50
G30652	DCB-D-50	Double-Ended	6.5-5	n/a	240	50
G30715	DCB-SV-50	Single-Ended	6.5-5	10, 6-4.5	240	50
G30714	DCB-DV-50	Double-Ended	6.5-5	10, 6-4.5	240	50

Minimum accessory channel 2.0 mm.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com

