

ShortShot*

SAEED HEMORRHOIDAL MULTI-BAND LIGATOR WITH TRIVIEW™

Unlike conventional devices, the ShortShot with TriView anoscope allows you to ligate up to three internal hemorrhoids during one procedure. With a trigger-based design for one-handed band deployment, the ShortShot resets quickly. And, since there is only a single insertion of the anoscope, patient discomfort is potentially reduced.





ShortShot*

SAEED HEMORRHOIDAL MULTI-BAND LIGATOR WITH TRIVIEW™

Used to ligate internal hemorrhoids facilitated by an anoscope.

Supplied non-sterile and is disposable - intended for single use only.

- Improved procedural efficiencies achieved through use of pre-loaded ShortShot Multi-Band Ligator
- Minimizes trauma with suction capture of hemorrhoids
- Superior positioning ability with TriView anoscope
- Safe, convenient and effective procedure
- Patient discomfort is potentially decreased with only a single insertion of an anoscope required

Order	Reference	Body	Diameter	Number
Number	Part Number	Length		of Bands
G34063	HMBL-4-TRI	5" / 12.7 cm	.6"/15 mm	4

For use with TriView Anoscope only. **Caution:** This product contains natural rubber latex, which may cause an allergic reaction. Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com Austria: +43 179567121, oe.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

www.cookmedical.com



Americas: EDI - www.cookmedical.com/edi.do

E-mail: orders@cookmedical.com

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335