
Manage the most challenging salivary stones with durable nitinol stone extractors.



Salivary Stone Extractors


Extracting salivary stones can be challenging. Cook Medical, a world leader in nitinol basket technologies, is delivering two salivary stone extractors that help make sialendoscopy more efficient and effective.

Nitinol Design


The shape memory and resilient qualities of our stone extractors' patented nitinol design allow expedient stone retrieval and minimize trauma to the ductal wall.

Features include:

- Kink resistance, which allows the baskets to withstand multiple deployments.
- Flexibility, which allows the expedient capture and release of stones.



NCircle®
Tipless Salivary Stone Extractor



NGage[®] Salivary Stone Extractor

NCircle®

TIPLESS SALIVARY STONE EXTRACTOR

Easily remove stones while reducing trauma to the ductal wall with the patented, tipless NCircle extractor.

Feature and Benefits

Tipless design

The tipless configuration allows you to position the basket directly against the mucosal lining, which helps minimize trauma, reduce the risk of ductal wall perforation and maintain a clear working field.

Increased flexibility

The flexible sheath allows continuous, full deflection of the endoscope during stone manipulation.



NGage® SALIVARY STONE EXTRACTOR

The NGage nitinol stone extractor is not just a grasper or stone basket; it's a combination of the two.

Feature and Benefits

Dual performance

NGage can hold stones like a basket and release them like a grasper, providing optimized control throughout the procedure.


Patented design

Its unique design allows you to engage, reposition, release and extract stones in various locations throughout the salivary duct.



Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Basket Diameter mm
NCircle Tipless Salivary Stone Extractor				
G23797	SSEB-1.5-115-10	1.5	115	10
NGage Salivary Stone Extractor				
G23798	SSEB-1.7-115-8	1.7	115	8

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.



Other products of interest

Along with our salivary stone extractors, Cook Medical offers salivary duct access tools, including wire guides, dilators and introducer sheaths that help you gain and maintain ductal access and make salivary stone extraction more efficient and effective. We are dedicated to the advancement of sialendoscopy and the prevention of unnecessary salivary gland removal.

Customer Service Centers

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fi.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com Switzerland - French: +41 448009609, fr.orders@cookmedical.com Switzerland - Italian: +41 448009609, it.orders@cookmedical.com Switzerland - German: +41 448009609, de.orders@cookmedical.com United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: cau.custserv @ cook medical.com

