



Expand your procedural options

with the first 8 mm x 8 cm balloon that can pass through a 5 Fr sheath.



Pliaform™ balloon texture



Starcast™ shoulders

Advance® 35LP

LOW-PROFILE PTA BALLOON DILATATION CATHETER

Exceptional rewrap

Starcast™ balloon shoulders facilitate easy sheath reinsertion so that you can repeat the treatment during the same procedure.

Reduced friction

The Pliaform™ balloon texture reduces friction by more than 50%* to help you cross lesions and retract the balloon after deflation.

Rapid inflation and deflation

The D-shaped lumen ensures rapid balloon inflation and deflation.

*Data on file.



www.cookmedical.com

ADVANCE 35LP LOW-PROFILE PTA DILATATION CATHETER

Order Number	Reference Part Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
80 cm Catheter Length				
G52325	PTA5-35-80-3-2.0	5	10	15
G52326	PTA5-35-80-3-4.0	5	10	15
G34348	PTA5-35-80-3-6.0	5	10	15
G52327	PTA5-35-80-3-8.0	5	10	15
G52328	PTA5-35-80-3-10.0	5	10	15
G52329	PTA5-35-80-3-12.0	5	10	15
G52330	PTA5-35-80-3-14.0	5	10	15
G34349	PTA5-35-80-3-17.0	5	10	15
G34350	PTA5-35-80-3-20.0	5	10	15
G34351	PTA5-35-80-4-1.5	5	10	15
G52331	PTA5-35-80-4-2.0	5	10	15
G34352	PTA5-35-80-4-2.5	5	10	15
G52332	PTA5-35-80-4-4.0	5	10	15
G34355	PTA5-35-80-4-6.0	5	10	15
G52333	PTA5-35-80-4-8.0	5	10	15
G52334	PTA5-35-80-4-10.0	5	10	15
G52335	PTA5-35-80-4-12.0	5	10	15
G52336	PTA5-35-80-4-14.0	5	10	15
G34353	PTA5-35-80-4-17.0	5	10	15
G34354	PTA5-35-80-4-20.0	5	10	15
G34356	PTA5-35-80-5-1.5	5	10	15
G52337	PTA5-35-80-5-2.0	5	10	15
G34357	PTA5-35-80-5-2.5	5	10	15
G52338	PTA5-35-80-5-4.0	5	10	15
G52339	PTA5-35-80-5-6.0	5	10	15
G52340	PTA5-35-80-5-8.0	5	10	15
G52341	PTA5-35-80-5-10.0	5	10	15
G52342	PTA5-35-80-5-12.0	5	10	15
G52343	PTA5-35-80-5-14.0	5	10	15
G34358	PTA5-35-80-5-17.0	5	10	15
G34359	PTA5-35-80-5-20.0	5	10	15
G34360	PTA5-35-80-6-1.5	5	10	15
G52344	PTA5-35-80-6-2.0	5	10	15
G34361	PTA5-35-80-6-2.5	5	10	15
G52345	PTA5-35-80-6-4.0	5	10	15
G52346	PTA5-35-80-6-6.0	5	10	15
G52347	PTA5-35-80-6-8.0	5	10	15
G52348	PTA5-35-80-6-10.0	5	10	15
G52349	PTA5-35-80-6-12.0	5	10	15
G52350	PTA5-35-80-6-14.0	5	10	15
G34362	PTA5-35-80-6-17.0	6	10	15
G34363	PTA5-35-80-6-20.0	6	10	15
G34364	PTA5-35-80-7-1.5	5	8	12

Order Number	Reference Part Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
80 cm Catheter Length				
G52351	PTA5-35-80-7-2.0	5	8	12
G34365	PTA5-35-80-7-2.5	5	8	12
G52352	PTA5-35-80-7-4.0	5	8	12
G52353	PTA5-35-80-7-6.0	5	8	12
G52354	PTA5-35-80-7-8.0	5	8	12
G52355	PTA5-35-80-7-10.0	6	8	12
G52356	PTA5-35-80-7-12.0	6	8	12
G52357	PTA5-35-80-7-14.0	6	8	12
G34366	PTA5-35-80-7-17.0	6	8	12
G34367	PTA5-35-80-7-20.0	6	8	12
G34368	PTA5-35-80-8-1.5	5	8	11
G52358	PTA5-35-80-8-2.0	5	8	11
G34369	PTA5-35-80-8-2.5	5	8	11
G52359	PTA5-35-80-8-4.0	5	8	11
G52360	PTA5-35-80-8-6.0	5	8	11
G52361	PTA5-35-80-8-8.0	5	8	11
G52362	PTA5-35-80-8-10.0	6	8	11
G52363	PTA5-35-80-8-12.0	6	8	11
G52364	PTA5-35-80-8-14.0	6	8	11
G34370	PTA5-35-80-8-17.0	6	8	11
G34371	PTA5-35-80-8-20.0	6	8	11
G34372	PTA5-35-80-9-2.0	6	8	11
G52365	PTA5-35-80-9-3.0	6	8	11
G52366	PTA5-35-80-9-4.0	6	8	11
G52367	PTA5-35-80-9-6.0	6	8	11
G52368	PTA5-35-80-9-8.0	6	8	11
G52369	PTA5-35-80-9-10.0	6	8	11
G52370	PTA5-35-80-9-12.0	6	8	11
G52371	PTA5-35-80-9-14.0	6	8	11
G34373	PTA5-35-80-10-2.0	6	8	11
G52372	PTA5-35-80-10-3.0	6	8	11
G52373	PTA5-35-80-10-4.0	6	8	11
G52374	PTA5-35-80-10-6.0	6	8	11
G52375	PTA5-35-80-10-8.0	6	8	11
G52376	PTA5-35-80-10-10.0	7	8	11
G52377	PTA5-35-80-10-12.0	7	8	11
G52378	PTA5-35-80-10-14.0	7	8	11
G34374	PTA5-35-80-12-2.0	7	5	8
G52379	PTA5-35-80-12-4.0	7	5	8
G52380	PTA5-35-80-12-8.0	7	5	8
G52381	PTA5-35-80-12-10.0	7	5	8
G52382	PTA5-35-80-12-12.0	7	5	8
G52383	PTA5-35-80-12-14.0	7	5	8

Order Number	Reference Part Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
135 cm Catheter Length				
G52245	PTA5-35-135-3-2.0	5	10	15
G52246	PTA5-35-135-3-4.0	5	10	15
G34375	PTA5-35-135-3-6.0	5	10	15
G52247	PTA5-35-135-3-8.0	5	10	15
G52248	PTA5-35-135-3-10.0	5	10	15
G52249	PTA5-35-135-3-12.0	5	10	15
G52250	PTA5-35-135-3-14.0	5	10	15
G34376	PTA5-35-135-3-17.0	5	10	15
G34377	PTA5-35-135-3-20.0	5	10	15
G34378	PTA5-35-135-4-1.5	5	10	15
G52251	PTA5-35-135-4-2.0	5	10	15
G34379	PTA5-35-135-4-2.5	5	10	15
G52252	PTA5-35-135-4-4.0	5	10	15
G35528	PTA5-35-135-4-6.0	5	10	15
G52253	PTA5-35-135-4-8.0	5	10	15
G52254	PTA5-35-135-4-10.0	5	10	15
G52255	PTA5-35-135-4-12.0	5	10	15
G52256	PTA5-35-135-4-14.0	5	10	15
G34380	PTA5-35-135-4-17.0	5	10	15
G35527	PTA5-35-135-4-20.0	5	10	15
G35529	PTA5-35-135-5-1.5	5	10	15
G52257	PTA5-35-135-5-2.0	5	10	15
G35530	PTA5-35-135-5-2.5	5	10	15
G52258	PTA5-35-135-5-4.0	5	10	15
G52259	PTA5-35-135-5-6.0	5	10	15
G52260	PTA5-35-135-5-8.0	5	10	15
G52261	PTA5-35-135-5-10.0	5	10	15
G52262	PTA5-35-135-5-12.0	5	10	15
G52263	PTA5-35-135-5-14.0	5	10	15
G35531	PTA5-35-135-5-17.0	5	10	15
G35532	PTA5-35-135-5-20.0	5	10	15
G35533	PTA5-35-135-6-1.5	5	10	15
G52264	PTA5-35-135-6-2.0	5	10	15
G35534	PTA5-35-135-6-2.5	5	10	15
G52265	PTA5-35-135-6-4.0	5	10	15
G52266	PTA5-35-135-6-6.0	5	10	15
G52267	PTA5-35-135-6-8.0	5	10	15
G52268	PTA5-35-135-6-10.0	5	10	15
G52269	PTA5-35-135-6-12.0	5	10	15
G52270	PTA5-35-135-6-14.0	5	10	15
G35535	PTA5-35-135-6-17.0	6	10	15
G35536	PTA5-35-135-6-20.0	6	10	15
G35537	PTA5-35-135-7-1.5	5	8	12
G52271	PTA5-35-135-7-2.0	5	8	12


Order Number	Reference Part Number	Compatible Sheath Fr	Nominal Pressure atm	Rated Burst Pressure atm
135 cm Catheter Length				
G35538	PTA5-35-135-7-2.5	5	8	12
G52272	PTA5-35-135-7-4.0	5	8	12
G52273	PTA5-35-135-7-6.0	5	8	12
G52274	PTA5-35-135-7-8.0	5	8	12
G52275	PTA5-35-135-7-10.0	6	8	12
G52276	PTA5-35-135-7-12.0	6	8	12
G52277	PTA5-35-135-7-14.0	6	8	12
G35539	PTA5-35-135-7-17.0	6	8	12
G35540	PTA5-35-135-7-20.0	6	8	12
G35541	PTA5-35-135-8-1.5	5	8	11
G52278	PTA5-35-135-8-2.0	5	8	11
G35542	PTA5-35-135-8-2.5	5	8	11
G52279	PTA5-35-135-8-4.0	5	8	11
G52280	PTA5-35-135-8-6.0	5	8	11
G52281	PTA5-35-135-8-8.0	5	8	11
G52282	PTA5-35-135-8-10.0	6	8	11
G52283	PTA5-35-135-8-12.0	6	8	11
G52284	PTA5-35-135-8-14.0	6	8	11
G35543	PTA5-35-135-8-17.0	6	8	11
G35544	PTA5-35-135-8-20.0	6	8	11
G35545	PTA5-35-135-9-2.0	6	8	11
G52285	PTA5-35-135-9-3.0	6	8	11
G52286	PTA5-35-135-9-4.0	6	8	11
G52287	PTA5-35-135-9-6.0	6	8	11
G52288	PTA5-35-135-9-8.0	6	8	11
G52289	PTA5-35-135-9-10.0	6	8	11
G52290	PTA5-35-135-9-12.0	6	8	11
G52291	PTA5-35-135-9-14.0	6	8	11
G35546	PTA5-35-135-10-2.0	6	8	11
G52292	PTA5-35-135-10-3.0	6	8	11
G52293	PTA5-35-135-10-4.0	6	8	11
G52294	PTA5-35-135-10-6.0	6	8	11
G52295	PTA5-35-135-10-8.0	6	8	11
G52296	PTA5-35-135-10-10.0	7	8	11
G52297	PTA5-35-135-10-12.0	7	8	11
G52298	PTA5-35-135-10-14.0	7	8	11
G35547	PTA5-35-135-12-2.0	7	5	8
G52299	PTA5-35-135-12-4.0	7	5	8
G52300	PTA5-35-135-12-8.0	7	5	8
G52301	PTA5-35-135-12-10.0	7	5	8
G52302	PTA5-35-135-12-12.0	7	5	8
G52303	PTA5-35-135-12-14.0	7	5	8

Some products or part numbers may not be available in all markets.
Contact your local Cook representative or Customer Service for details.


The following products are compatible with Advance 35LP.




The **Zilver® PTX® Drug-Eluting Peripheral Stent** is the world's only drug-eluting stent for the superficial femoral artery.




The **Amplatz Support Wire Guide** is available in three stiffness options and is designed for atraumatic navigation and support in tortuous anatomy.



The **Flexor® Ansel Guiding Sheath with High-Flex Dilator** features our most advanced technologies for enhanced pushability and trackability and is available in 4, 5, 6, 7, 8, 10, and 12 Fr.



The **Slip-Cath® Angiographic Catheter** features a lubricious coating that can help you to smoothly navigate tortuous anatomy.



The **CXI® Support Catheter** provides the strength, visibility, lubricity, and low profile that you need to help you cross lesions.

Pressure (bar) (atm)	Balloon Diameter (mm)								
	3	4	5	6	7	8	9	10	12
4	2.64	3.49	4.28	5.26	6.31	7.17	8.08	9.11	11.62
5	2.69	3.58	4.40	5.39	6.52	7.42	8.37	9.40	12.00
6	2.75	3.67	4.54	5.55	6.72	7.66	8.62	9.65	12.32
7	2.82	3.77	4.67	5.68	6.88	7.85	8.84	9.84	12.58
8	2.87	3.86	4.79	5.80	7.00	8.00	9.00	10.00	12.84
9	2.95	3.93	4.90	5.90	7.12	8.15	9.15	10.16	
10	3.00	4.00	5.00	6.00	7.23	8.27	9.33	10.29	
11	3.05	4.06	5.09	6.09	7.34	8.41	9.49	10.47	
12	3.10	4.12	5.19	6.17	7.46				
13	3.15	4.18	5.28	6.26					
14	3.20	4.24	5.37	6.34					
15	3.25	4.31	5.47	6.44					

= Nominal Balloon Pressure
 = Rated Burst Pressure (RBP)

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com

