

Occlude a wide range
of vessel sizes with
coil lengths from
3-20 cm and diameters
from 2-20 mm.



Now
available
in more
sizes.

Nester[®]
EMBOLIZATION COIL

**Expand your options with
a large number of coil profiles**

Use with a complementary line of
embolization products to support
your procedure from start to finish.

Excellent packing coil

Soft platinum coils pack densely with
soft fiber for excellent occlusion.

Use alone or in combination

Use in conjunction with other
pushable coils, vascular plugs and
detachable coil systems.



www.cookmedical.com

NOW AVAILABLE IN A BROADER SELECTION OF MICROCOILS

Rows shaded in blue are new additions to the product line.

	Order Number	Reference Part Number	Recommended Catheter ID/ End Hole Diameter inch	Extended Embolus Length cm	Coiled Embolus Diameter mm	Approximate Number of Loops
.018 Inch Microcoils™						
3cm Extended Embolus	G52731	MWCE-18-3-2-NESTER	.018	3	2	4.8
	G52732	MWCE-18-3-3-NESTER	.018	3	3	3.2
5cm Extended Embolus	G52733	MWCE-18-5-2-NESTER	.018	5	2	8.0
	G52734	MWCE-18-5-3-NESTER	.018	5	3	5.3
	G23520	MWCE-18-5-5-NESTER	.018	5	5	3.2
7cm Extended Embolus	G52735	MWCE-18-7-2-NESTER	.018	7	2	11.1
	G52736	MWCE-18-7-3-NESTER	.018	7	3	7.4
	G52737	MWCE-18-7-4-NESTER	.018	7	4	5.6
	G23521	MWCE-18-7-5-NESTER	.018	7	5	4.5
	G52738	MWCE-18-7-6-NESTER	.018	7	6	3.7
	G52739	MWCE-18-7-8-NESTER	.018	7	8	2.8
	G52741	MWCE-18-7-10-NESTER	.018	7	10	2.2
	G26987	MWCE-18-14-3-NESTER	.018	14	3	14.9
14cm Extended Embolus	G26988	MWCE-18-14-4-NESTER	.018	14	4	11.1
	G23522	MWCE-18-14-5-NESTER	.018	14	5	8.9
	G26989	MWCE-18-14-6-NESTER	.018	14	6	7.4
	G26990	MWCE-18-14-8-NESTER	.018	14	8	5.6
	G26986	MWCE-18-14-10-NESTER	.018	14	10	4.5

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

NOW AVAILABLE IN A BROADER SELECTION OF MACROCOILS

Rows shaded in blue are new additions to the product line.

	Order Number	Reference Part Number	Recommended Catheter ID/ End Hole Diameter inch	Extended Embolus Length cm	Coiled Embolus Diameter mm	Approximate Number of Loops
.035 Inch Coils						
7cm Extended Embolus	G23523	MWCE-35-7-3-NESTER	.035	7	3	7.4
	G52742	MWCE-35-7-4-NESTER	.035	7	4	5.6
	G23524	MWCE-35-7-5-NESTER	.035	7	5	4.5
	G52743	MWCE-35-7-6-NESTER	.035	7	6	3.7
	G52744	MWCE-35-7-8-NESTER	.035	7	8	2.8
	G52745	MWCE-35-7-10-NESTER	.035	7	10	2.2
	G52746	MWCE-35-7-12-NESTER	.035	7	12	1.9
	G52747	MWCE-35-7-14-NESTER	.035	7	14	1.6
	G52748	MWCE-35-7-16-NESTER	.035	7	16	1.4
	G52749	MWCE-35-7-18-NESTER	.035	7	18	1.2
	G52750	MWCE-35-7-20-NESTER	.035	7	20	1.1
14cm Extended Embolus	G23525	MWCE-35-14-3-NESTER	.035	14	3	14.9
	G26993	MWCE-35-14-4-NESTER	.035	14	4	11.1
	G23526	MWCE-35-14-5-NESTER	.035	14	5	8.9
	G26994	MWCE-35-14-6-NESTER	.035	14	6	7.4
	G26995	MWCE-35-14-8-NESTER	.035	14	8	5.6
	G26991	MWCE-35-14-10-NESTER	.035	14	10	4.5
	G26992	MWCE-35-14-12-NESTER	.035	14	12	3.7
	G52751	MWCE-35-14-14-NESTER	.035	14	14	3.2
	G52752	MWCE-35-14-16-NESTER	.035	14	16	2.8
	G52753	MWCE-35-14-18-NESTER	.035	14	18	2.5
	G52754	MWCE-35-14-20-NESTER	.035	14	20	2.2
20cm Extended Embolus	G23527	MWCE-35-20-5-NESTER	.035	20	5	12.7
	G23528	MWCE-35-20-6-NESTER	.035	20	6	10.6
	G23529	MWCE-35-20-8-NESTER	.035	20	8	8
	G23530	MWCE-35-20-10-NESTER	.035	20	10	6.4
	G23531	MWCE-35-20-12-NESTER	.035	20	12	5.3
	G23532	MWCE-35-20-14-NESTER	.035	20	14	4.5
	G23533	MWCE-35-20-16-NESTER	.035	20	16	4.0
	G52755	MWCE-35-20-18-NESTER	.035	20	18	3.5
	G52756	MWCE-35-20-20-NESTER	.035	20	20	3.2

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Customer Service

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 972519996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: cau.custserv@cookmedical.com

