




Slip-Cath[®]

ANGIOGRAPHIC CATHETERS



Trusted performance in all
of your procedures, from the
everyday to the extreme.



Slip-Cath[®]


ANGIOGRAPHIC CATHETERS

Superior ease of insertion

Slip-Cath's AQ[®] hydrophilic coating is 5 times more lubricious than noncoated catheters. The tapered tip slides easily into selected vessels with minimal force and less trauma.

Superior trackability

Slip-Cath's innovative design delivers superior torque control to the catheter tip, ensuring maximum trackability through tortuous anatomy.




Superior visualization

Cook's Beacon[®] enhanced radiopaque tip technology ensures excellent visibility during catheter manipulation.

Superior kink resistance

Nylon catheter construction with stainless steel braiding from hub to tip provides maximum control and pushability, minimizing kinking in tight curves.



Catheter Flow Rates


Flow rates shown are determined in vitro using the Cook Pressure Injector and 100% Oxilan® 350 at room temperature (25°C). Increasing the temperature of the contrast medium (i.e., to 37°C) or using a less viscous contrast medium will result in an increase in flow rate. The viscosity of Oxilan 350 at room temperature (25°C) is 16.3 cP vs Renografin-76, which at room temperature (25°C) is 18.5 cP. Maximum injection pressure is 1,200 psi (82 kg/cm²). Since a number of variables can produce different flow rates, the following information should be used only as a guide.

SCBR4.0-38				
French Size	Length	MAXIMUM FLOW RATE	MAXIMUM RECOMMENDED PRESSURE	
		cc/sec		
		100% Oxilan 350	kg/cm ²	psi
4.0	40	25	82	1,200
4.0	65	21	82	1,200
4.0	80	19	82	1,200
4.0	100	15	82	1,200




SCBR5.0-38				
French Size	Length	MAXIMUM FLOW RATE	MAXIMUM RECOMMENDED PRESSURE	
		cc/sec		
		100% Oxilan 350	kg/cm ²	psi
5.0	40	25	82	1,200
5.0	65	21	82	1,200
5.0	80	19	82	1,200
5.0	90	17	82	1,200
5.0	100	15	82	1,200
5.0	125	11	82	1,200

Oxilan is a registered trademark of Guerbet LLC, Bloomington, Indiana.




French/Diameters



Cerebral





Visceral




Lower Extremity





Contralateral Approach



Uterine



Multipurpose






CEREBRAL


Global Product Number	Order Number	French Size	Accepts Wire Guide Diameter inch	Length cm
G31016	SCBR4.0-38-100-P-NS-H1	4.0	.038	100
G11991	SCBR5.0-38-100-P-NS-H1	5.0	.038	100
G31327	SCBR4.0-38-100-P-NS-JB1	4.0	.038	100
G11997	SCBR5.0-38-100-P-NS-JB1	5.0	.038	100
G31328	SCBR4.0-38-100-P-NS-JB2	4.0	.038	100
G11998	SCBR5.0-38-100-P-NS-JB2	5.0	.038	100
G31325	SCBR4.0-38-100-P-NS-SIM1	4.0	.038	100
G12000	SCBR5.0-38-100-P-NS-SIM1	5.0	.038	100
G31017	SCBR4.0-38-100-P-NS-SIM2	4.0	.038	100
G11992	SCBR5.0-38-100-P-NS-SIM2	5.0	.038	100
G13163	SCBR5.0-38-100-P-NS-VTK	5.0	.038	100
G12082	SCBR5.0-38-125-P-NS-VTK	5.0	.038	125
G12168	SCBR5.0-38-125-P-NS-VERT	5.0	.038	125
G31330	SCBR4.0-38-100-P-NS-VERT	4.0	.038	100
G12007	SCBR5.0-38-100-P-NS-VERT	5.0	.038	100
G12043	SCR3.0B-NT-120-P-NS-0-RT	3.0	.025	120
G12044	SCR3.0B-NT-135-P-NS-0-RT	3.0	.025	135
G12045	SCR3.0B-NT-150-P-NS-0-RT	3.0	.025	150

Order Number Key

SCBR4.0-38-100-P-NS-H1


- Tip Shape
- 100=Catheter Length (cm)
- 38=Wire Guide Diameter (.0XX inch)
- 4.0=French Size






VISCERAL

Global Product Number	Order Number	French Size	Accepts Wire Guide Diameter inch	Length cm
G12003	SCBR5.0-38-65-P-NS-C1	5.0	.038	65
G31021	SCBR4.0-38-65-P-NS-C2	4.0	.038	65
G11994	SCBR5.0-38-65-P-NS-C2	5.0	.038	65
G12006	SCBR5.0-38-80-P-NS-C2	5.0	.038	80
G31022	SCBR4.0-38-100-P-NS-C2	4.0	.038	100
G12008	SCBR5.0-38-100-P-NS-C2	5.0	.038	100
G31004	SCBR4.0-38-100-P-NS-DAV	4.0	.038	100
G11996	SCBR5.0-38-100-P-NS-DAV	5.0	.038	100
G31018	SCBR4.0-38-40-P-NS-KMP	4.0	.038	40
G11993	SCBR5.0-38-40-P-NS-KMP	5.0	.038	40
G31019	SCBR4.0-38-65-P-NS-KMP	4.0	.038	65
G12156	SCBR5.0-38-65-P-NS-KMP	5.0	.038	65
G31020	SCBR4.0-38-100-P-NS-KMP	4.0	.038	100
G31023	SCBR4.0-38-100-P-NS-MPA	4.0	.038	100
G12002	SCBR5.0-38-100-P-NS-MPA	5.0	.038	100
G36214	SCBR5.0-38-80-P-NS-RH	5.0	.038	80
G11995	SCBR5.0-38-65-P-NS-RIM	5.0	.038	65
G36319	SCBR4.0-38-65-P-NS-TC	4.0	.038	65
G36320	SCBR5.0-38-65-P-NS-TC	5.0	.038	65
G36321	SCBR4.0-38-65-P-NS-TC-BNK	4.0	.038	65
G36322	SCBR5.0-38-65-P-NS-TC-BNK	5.0	.038	65
G12179	SCBR5.0-38-65-P-NS-TEGT	5.0	.038	65
G31330	SCBR4.0-38-100-P-NS-VERT	4.0	.038	100
G12007	SCBR5.0-38-100-P-NS-VERT	5.0	.038	100
G12168	SCBR5.0-38-125-P-NS-VERT	5.0	.038	125
G13431	SCBR5.0-38-80-P-NS-VS	5.0	.038	80
G13432	SCBR5.0-38-80-P-NS-VS1	5.0	.038	80
G13433	SCBR5.0-38-80-P-NS-VS2	5.0	.038	80
G31335	SCBR4.0-38-65-P-NS-0	4.0	.038	65
G12009	SCBR5.0-38-65-P-NS-0	5.0	.038	65
G31336	SCBR4.0-38-100-P-NS-0	4.0	.038	100
G12005	SCBR5.0-38-100-P-NS-0	5.0	.038	100
G12046	SCR3.0B-NT-80-P-NS-0-RT	3.0	.025	80
G12042	SCR3.0B-NT-100-P-NS-0-RT	3.0	.025	100
G12043	SCR3.0B-NT-120-P-NS-0-RT	3.0	.025	120
G12044	SCR3.0B-NT-135-P-NS-0-RT	3.0	.025	135
G12045	SCR3.0B-NT-150-P-NS-0-RT	3.0	.025	150






LOWER EXTREMITY


Global Product Number	Order Number	French Size	Accepts Wire Guide Diameter inch	Length cm
G31335	SCBR4.0-38-65-P-NS-0	4.0	.038	65
G12009	SCBR5.0-38-65-P-NS-0	5.0	.038	65
G31336	SCBR4.0-38-100-P-NS-0	4.0	.038	100
G12005	SCBR5.0-38-100-P-NS-0	5.0	.038	100
G36391	SCBR4.0-35-150-P-NS-0	4.0	.035	150
G12046	SCR3.0B-NT-80-P-NS-0-RT	3.0	.025	80
G12042	SCR3.0B-NT-100-P-NS-0-RT	3.0	.025	100
G12043	SCR3.0B-NT-120-P-NS-0-RT	3.0	.025	120
G12044	SCR3.0B-NT-135-P-NS-0-RT	3.0	.025	135
G12045	SCR3.0B-NT-150-P-NS-0-RT	3.0	.025	150



CONTRALATERAL APPROACH


Global Product Number	Order Number	French Size	Accepts Wire Guide Diameter inch	Length cm
G31021	SCBR4.0-38-65-P-NS-C2	4.0	.038	65
G11994	SCBR5.0-38-65-P-NS-C2	5.0	.038	65
G12006	SCBR5.0-38-80-P-NS-C2	5.0	.038	80
G31022	SCBR4.0-38-100-P-NS-C2	4.0	.038	100
G12008	SCBR5.0-38-100-P-NS-C2	5.0	.038	100
G13431	SCBR5.0-38-80-P-NS-VS	5.0	.038	80
G13432	SCBR5.0-38-80-P-NS-VS1	5.0	.038	80
G13433	SCBR5.0-38-80-P-NS-VS2	5.0	.038	80
G11995	SCBR5.0-38-65-P-NS-RIM	5.0	.038	65
G36319	SCBR4.0-38-65-P-NS-TC	4.0	.038	65
G36320	SCBR5.0-38-65-P-NS-TC	5.0	.038	65
G36321	SCBR4.0-38-65-P-NS-TC-BNK	4.0	.038	65
G36322	SCBR5.0-38-65-P-NS-TC-BNK	5.0	.038	65





UTERINE

Global Product Number	Order Number	French Size	Accepts Wire Guide Diameter inch	Length cm
G12006	SCBR5.0-38-80-P-NS-C2	5.0	.038	80
G31022	SCBR4.0-38-100-P-NS-C2	4.0	.038	100
G12008	SCBR5.0-38-100-P-NS-C2	5.0	.038	100
G31234	SCBR5.0-35-90-P-NS-RUC	5.0	.035	90
G36319	SCBR4.0-38-65-P-NS-TC	4.0	.038	65
G36320	SCBR5.0-38-65-P-NS-TC	5.0	.038	65
G36321	SCBR4.0-38-65-P-NS-TC-BNK	4.0	.038	65
G36322	SCBR5.0-38-65-P-NS-TC-BNK	5.0	.038	65
G12046	SCR3.0B-NT-80-P-NS-0-RT	3.0	.025	80
G12042	SCR3.0B-NT-100-P-NS-0-RT	3.0	.025	100
G12043	SCR3.0B-NT-120-P-NS-0-RT	3.0	.025	120





MULTIPURPOSE

Global Product Number	Order Number	French Size	Accepts Wire Guide Diameter inch	Length cm
G31004	SCBR4.0-38-100-P-NS-DAV	4.0	.038	100
G11996	SCBR5.0-38-100-P-NS-DAV	5.0	.038	100
G31018	SCBR4.0-38-40-P-NS-KMP	4.0	.038	40
G11993	SCBR5.0-38-40-P-NS-KMP	5.0	.038	40
G31019	SCBR4.0-38-65-P-NS-KMP	4.0	.038	65
G12156	SCBR5.0-38-65-P-NS-KMP	5.0	.038	65
G31020	SCBR4.0-38-100-P-NS-KMP	4.0	.038	100
G31023	SCBR4.0-38-100-P-NS-MPA	4.0	.038	100
G12002	SCBR5.0-38-100-P-NS-MPA	5.0	.038	100
G12179	SCBR5.0-38-65-P-NS-TEGT	5.0	.038	65
G31330	SCBR4.0-38-100-P-NS-VERT	4.0	.038	100
G12007	SCBR5.0-38-100-P-NS-VERT	5.0	.038	100
G12168	SCBR5.0-38-125-P-NS-VERT	5.0	.038	125



DAV



KMP



MPA



TEGT



VERT



www.cookmedical.com

COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.

Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA

Phone: 905.640.7110, Toll Free: 800.668.0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.

Brisbane Technology Park, 12 Electronics Street, Eight Mile Plains

Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

WILLIAM COOK EUROPE ApS

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2008

PI-BM-SCDS-EN-200811

AORTIC INTERVENTION

CARDIOLOGY

CRITICAL CARE

ENDOSCOPY

PERIPHERAL INTERVENTION

SURGERY

UROLOGY

WOMEN'S HEALTH