




ORDERING INFORMATION



Used for the administration of nutrients, chemotherapeutic agents and other drugs for therapy, and for blood sampling, blood delivery and venous pressure monitoring.

- Cook Spectrum catheters are impregnated with the antibiotics minocycline and rifampin, which help to provide protection against CRBSIs.
- The Peel-Away sheath allows percutaneous introduction of the silicone catheter.
- 21 gage venous bevel echogenic needle.

Supplied sterile in peel-open packages. Intended for one-time use.

Global Product Number	Order Number	Catheter Fr	Catheter Length cm	Sheath Length cm	Wire Guide Diameter inch	Wire Length cm				
Single Lumen										
G51196	HASC-4.0-ABRM	4	12	9	.018	40				
G51197	HASC-5.0-ABRM	5	25	10	.018	40				
Double Lumen										
G51198	HASDC-7.0-LSC-ABRM	7	25	9	.018	40				
G51199	HASDC-7.0-RSC-ABRM	7	20	9	.018	40				
G51201	HASDC-9.0-ABRM	9	25	15.5	.018	60				
Triple Lumen										
G51200	HASTC-8.0-ABRM	8	25	15.5	.018	60				

Set Components: Catheter, wire, Peel-Away sheath, 21 gage needle, disposable syringe, disposable scalpel, needle holder, measuring tape and CLC2000 caps.

Lumen and Flow Rate Information

	Fr	Lumen Volume cc		Inner Diameter inch		Lumen Shape		Initial Length cm	Max Flow Rate cc/min	
Single Lumen	4	.1		.026		circle		12	16	
	5	.2		.036		circle		25	25	
		Small	Large	Small	Large	Small	Large		Small	Large
Double Lumen	7 LSC*	.05	.06	.019 x .045	.025 x .050	D	D	25	13	30
	7 RSC*	.05	.06	.019 x .045	.025 x .051	D	D	20	13	30
	9	1	1.1	.025 x .065	.036 x .074	D	D	25	44	108
		Small (2)	Large	Small (2)	Large	Small (2)	Large		Small (2)	Large
Triple Lumen	8	.9	1.1	.029	.043	circle	circle	25	15	50

*LSC=Left subclavian

*RSC=Right subclavian



www.cookmedical.com

COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.

Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA

Phone: 905.640.7110, Toll Free: 800.668.0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.

Brisbane Technology Park, 12 Electronics Street, Eight Mile Plains Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 38 41 11 88

WILLIAM COOK EUROPE ApS

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2008

PI-BM-SHCDS-EN-200810

AORTIC INTERVENTION | CARDIOLOGY | CRITICAL CARE | ENDOSCOPY | PERIPHERAL INTERVENTION | SURGERY | UROLOGY | WOMEN'S HEALTH