




Wire Guides

CATALOG



Key for Anatomy Blocks



Contents

HYDROPHILIC WIRE GUIDES

Roadrunner® Hydrophilic Wire Guide	5
HiWire® Hydrophilic Wire Guide	6

CLASSIC WIRE GUIDES

Fixed-Core Wire Guide	9
Bentson Fixed-Core Wire Guide	11
Newton Fixed-Core Wire Guide	12
Rosen Fixed-Core Wire Guide	13
Double Flexible-Tipped Fixed-Core Wire Guide	14
Movable Core Wire Guide	15


SPECIALTY WIRE GUIDES


Approach® CTO Microwire Guide	17
Approach® Hydro ST Microwire Guide	18
Approach® Pro-LT Microwire Guide	18
Shorty® Interventional Wire Guide	19
Coons Interventional Wire Guide	19
Roadrunner® Extra-Support Mandril Wire Guide	20
Graduate™ Measuring Wire Guide	20
McNamara Renal Exchange Wire Guide	21
Lunderquist-Ring Torque Interventional Wire Guide	21
Reuter Tip-Deflecting Wire Guide	22
Cope Mandril Wire Guide	23
Torq-Flex® Mandril Wire Guide	23

SUPPORT WIRE GUIDES

Amplatz Support Wire Guide	25
Lunderquist® Support Wire Guide	26

Reference Part Number Key





Hydrophilic Wire Guides


Hydrophilic wire guides offer physicians increased trackability and pushability in both simple and complex cases. These nitinol-core wires feature technologies that provide smooth navigation and maximize control and visibility—even in tortuous anatomy.

Roadrunner®

HYDROPHILIC WIRE GUIDE

Used for complex diagnostic and interventional procedures.

- **AQ® hydrophilic coating** reduces friction and provides optimal trackability and pushability
- **Soft, radiopaque platinum** tips allow atraumatic vessel selection
- **Specially tapered nitinol mandril** provides superb directional control and 1:1 torque response



Order Number	Reference Part Number
Nimble®	
Short Taper/Flexible Shaft	
G09607	RFSPC-35-145
G09608	RFSPC-35-180
G09609	RFSPC-35-260
G09610	RFSPC-38-145
Nimble Floppy	
Standard Taper/Flexible Shaft	
G07914	RFPC-35-145
G07937	RFPC-35-180
The Firm™	
Standard Taper/Stiff Shaft	
G09534	RPC-35-80
G06979	RPC-35-145
G07285	RPC-35-145-0
G07518	RPC-35-180
G09105	RPC-35-260
G06980	RPC-38-145


Order Number	Reference Part Number
The Firm LT	
Long Taper/Stiff Shaft	
G06866	RLPC-35-145
G07287	RLPC-35-145-0
G07516	RLPC-35-180
Nimble – Connors Modification	
Short Taper/Flexible Shaft	
G12150	RFSPC-35-145-CONNORS
The Firm – Connors Modification	
Standard Taper/Stiff Shaft	
G11666	RPC-35-145-CONNORS
G11667	RPC-35-260-2-CONNORS

HiWire®

HYDROPHILIC WIRE GUIDE

Used to facilitate the placement of devices during diagnostic and interventional procedures.

- Enhanced coating provides **long-lasting lubricity**
- Optimum shaft stiffness and tip flexibility, combined with 1:1 torque response, **provide total control**
- Advanced radiopaque materials enhance the visibility of HiWire hydrophilic wire guides




Order Number	Reference Part Number
Standard Shaft	
Straight Tip	
G36276	HW-18-180
G36334	HW-25-150
G36283	HW-35-150
G36324	HW-35-180
G36325	HW-35-260
G36287	HW-38-150
G36335	HW-38-180
Standard Shaft	
Angled Tip	
G48670	HWA-18-80
G36275	HWA-18-150
G36277	HWA-18-180
G36278	HWA-18-260
G36279	HWA-25-150
G36280	HWA-25-180
G36281	HWA-25-260
G48671	HWA-35-80
G36282	HWA-35-150
G36284	HWA-35-180
G36300	HWA-35-260
G48672	HWA-38-80
G36286	HWA-38-150
G36288	HWA-38-180
G36301	HWA-38-260

Order Number	Reference Part Number
Stiff Shaft	
Straight Tip	
G36291	HWS-35-150
G36326	HWS-35-180
G36327	HWS-35-260
G36295	HWS-38-150
G36336	HWS-38-180
Stiff Shaft	
Angled Tip	
G48673	HWAS-35-80
G36290	HWAS-35-150
G36292	HWAS-35-180
G36293	HWAS-35-260
G48674	HWAS-38-80
G36294	HWAS-38-150
G36296	HWAS-38-180
G36297	HWAS-38-260

HiWire®

HYDROPHILIC WIRE GUIDE



Get there with versatile hydrophilic technology.

ENHANCED COATING


New and improved coating provides long-lasting lubricity to facilitate diagnostic and interventional procedures.

TOTAL CONTROL

Optimum shaft stiffness and tip flexibility, combined with 1:1 torque response, provide maximum directional control.

GREATER VISIBILITY

Advanced materials enhance the visibility of HiWire hydrophilic wire guides.



Classic Wire Guides


Through partnerships with pioneering physicians, Cook Medical developed many wire guide technologies that are now industry standards. These fixed- and movable-core, PTFE-coated, stainless steel wire guides are classics that interventionalists have been relying on for years. Available in a variety of options, the classics are ideal for everyday procedures.

Fixed-Core

WIRE GUIDE

Used for vascular and nonvascular applications.

- Tapered fixed-core construction provides **gradual transition to a flexible tip**
- Available with **Safe-T-J® curved tip** to minimize vessel trauma
- PTFE-coated stainless steel




Order Number	Reference Part Number
Standard Straight	
PTFE-Coated Stainless Steel	
G00593	TSF-18-50
G00587	TSF-18-145
G00590	TSF-18-260
G00609	TSF-21-50
G00602	TSF-21-145
G00605	TSF-21-260
G00624	TSF-25-50
G00626	TSF-25-80
G00615	TSF-25-125
G00617	TSF-25-145
G00620	TSF-25-260
G00638	TSF-32-145
G00661	TSF-35-50
G00664	TSF-35-80
G00650	TSF-35-145
G00652	TSF-35-180
G00655	TSF-35-260
G00670	TSF-38-145
G00674	TSF-38-260
Standard Straight	
PTFE-Coated Stainless Steel with Heparin Coating	
G00588	TSF-18-145-BH
G00603	TSF-21-145-BH
G00618	TSF-25-145-BH
G00639	TSF-32-145-BH
G00662	TSF-35-50-BH
G00651	TSF-35-145-BH
G00656	TSF-35-260-BH
G00671	TSF-38-145-BH

Order Number	Reference Part Number
Heavy Duty Straight	
PTFE-Coated Stainless Steel	
G02109	THSF-25-180
G01057	THSF-35-145
G01518	THSF-38-80
G00426	THSF-38-145
Heavy Duty Straight	
PTFE-Coated Stainless Steel with Short Taper	
G06431	THSF-35-180-NY
Heavy Duty Straight	
PTFE-Coated Stainless Steel with Heparin Coating	
G00422	THSF-35-145-BH
G00427	THSF-38-145-BH
Standard Safe-T-J	
PTFE-Coated Stainless Steel	
G00769	TSCF-18-50-3
G00436	TSCF-18-145-3
G00456	TSCF-21-50-3
G00452	TSCF-21-145-3
G28824	TSCF-21-160-3
G01793	TSCF-21-180-3
G00474	TSCF-25-50-3
G00466	TSCF-25-145-3
G01602	TSCF-25-180-3
G00471	TSCF-25-260-3
G00488	TSCF-32-145-3
G02335	TSCF-32-180-3
G00492	TSCF-32-260-3
G00524	TSCF-35-50-3
G00529	TSCF-35-80-3
G00527	TSCF-35-80-15

Fixed-Core (continued)

WIRE GUIDE

Order Number	Reference Part Number
Standard Safe-T-J (continued)	
PTFE-Coated Stainless Steel	
G00507	TSCF-35-145-1.5
G00511	TSCF-35-145-3
G00513	TSCF-35-145-7.5
G00509	TSCF-35-145-15
G00902	TSCF-35-180-3
G02481	TSCF-35-180-15
G00515	TSCF-35-200-3
G00517	TSCF-35-260-3
G00516	TSCF-35-260-15
G00550	TSCF-38-50-3
G01285	TSCF-38-60-3
G00552	TSCF-38-80-3
G00533	TSCF-38-100-3
G00541	TSCF-38-145-3
G00543	TSCF-38-145-7.5
G00539	TSCF-38-145-15
G00791	TSCF-35-260-1.5
G00546	TSCF-38-260-3
G00548	TSCF-38-260-7.5
Standard Safe-T-J	
PTFE-Coated Stainless Steel with Heparin Coating	
G00437	TSCF-18-145-3-BH
G00453	TSCF-21-145-3-BH
G01363	TSCF-21-180-3-BH
G00467	TSCF-25-145-3-BH
G02848	TSCF-25-180-3-BH
G00472	TSCF-25-260-3-BH
G00489	TSCF-32-145-3-BH
G00525	TSCF-35-50-3-BH
G00530	TSCF-35-80-3-BH
G00907	TSCF-35-145-1-BH
G00508	TSCF-35-145-1.5-BH
G00512	TSCF-35-145-3-BH
G00514	TSCF-35-145-7.5-BH
G00510	TSCF-35-145-15-BH
G00916	TSCF-35-180-3-BH
G00846	TSCF-35-200-3-BH
G00518	TSCF-35-260-3-BH
G00553	TSCF-38-80-3-BH



Order Number	Reference Part Number
G00542	TSCF-38-145-3-BH
G00540	TSCF-38-145-15-BH
G00547	TSCF-38-260-3-BH

Heavy Duty Safe-T-J

PTFE-Coated Stainless Steel

G01654	THSCF-18-145-3
G00869	THSCF-25-145-3
G02165	THSCF-25-180-3
G02384	THSCF-25-260-3
G01512	THSCF-35-80-1.5
G04069	THSCF-35-80-3
G01725	THSCF-35-80-15
G00412	THSCF-35-145-3
G00411	THSCF-35-145-15
G01151	THSCF-35-180-15
G01197	THSCF-35-260-3
G02800	THSCF-35-260-7.5
G01827	THSCF-38-80-3
G07763	THSCF-38-80-7.5
G01886	THSCF-38-80-15
G00414	THSCF-38-145-3
G01149	THSCF-38-145-15
G01348	THSCF-38-180-3
G01139	THSCF-38-180-15
G01196	THSCF-38-260-3
G02989	THSCF-38-260-7.5

Heavy Duty Safe-T-J

PTFE-Coated Stainless Steel with Heparin Coating

G00964	THSCF-35-145-3-BH
--------	-------------------


Not all part numbers shown on this product sheet may be approved for sale in all jurisdictions.
Consult with your local Cook representative or customer service center for details.

Bentson

FIXED-CORE WIRE GUIDE

Used for vascular and nonvascular applications.

- Atraumatic tip with an **extra-floppy 6 cm distal portion**
- **Distal flexibility** for negotiating tortuous areas of the anatomy
- **Large-diameter core** provides additional body for catheter control
- PTFE-coated stainless steel



Order Number	Reference Part Number
--------------	-----------------------

Standard

PTFE-Coated Stainless Steel

G04194	TSFB-18-180
G00856	TSFB-25-145
G00688	TSFB-32-145
G01315	TSFB-35-80
G01282	TSFB-35-100
G00691	TSFB-35-145
G01290	TSFB-35-180
G01095	TSFB-35-200
G01063	TSFB-35-260
G00693	TSFB-38-145

Standard

PTFE-Coated Stainless Steel with Heparin Coating

G01750	TSFB-35-80-BH
G00692	TSFB-35-145-BH
G01431	TSFB-35-180-BH
G00793	TSFB-38-145-BH

Order Number	Reference Part Number
--------------	-----------------------

Plus

PTFE-Coated Stainless Steel

G09735	TSFBP-35-145
G09782	TSFBP-35-180
G10126	TSFBP-35-260

Plus

PTFE-Coated Stainless Steel with Heparin Coating


G09783	TSFBP-35-145-BH
--------	-----------------

Newton

FIXED-CORE WIRE GUIDE

Used for vascular and nonvascular applications.

- **Atraumatic, flexible tip** for negotiating tortuous anatomy
- Available with **10 cm or 15 cm long tapered core**
- PTFE-coated stainless steel



Order Number	Reference Part Number
Straight LT – 10 cm Tapered Core PTFE-Coated Stainless Steel	
G00701	TSFNA-35-145
G00794	TSFNA-35-260
G00708	TSFNA-38-145
Straight LT – 10 cm Tapered Core PTFE-Coated Stainless Steel with Heparin Coating	
G00702	TSFNA-35-145-BH
G00709	TSFNA-38-145-BH
Straight LLT – 15 cm Tapered Core PTFE-Coated Stainless Steel	
G00715	TSFNB-35-145
G00719	TSFNB-38-145
Straight LLT – 15 cm Tapered Core PTFE-Coated Stainless Steel with Heparin Coating	
G00716	TSFNB-35-145-BH
G00720	TSFNB-38-145-BH
Curved LT – 10 cm Tapered Core PTFE-Coated Stainless Steel	
G00557	TSCFNA-35-145-1.5
G00561	TSCFNA-35-145-3
G00559	TSCFNA-35-145-15


Order Number	Reference Part Number
Curved LT – 10 cm Tapered Core PTFE-Coated Stainless Steel with Heparin Coating	
G00562	TSCFNA-35-145-3-BH
G00563	TSCFNA-35-145-7.5-BH
G00560	TSCFNA-35-145-15-BH
G00568	TSCFNA-38-145-3-BH
Curved LLT – 15 cm Tapered Core PTFE-Coated Stainless Steel	
G00574	TSCFNB-35-145-3
G00576	TSCFNB-35-145-7.5
G00572	TSCFNB-35-145-15
G00580	TSCFNB-38-145-3
G00578	TSCFNB-38-145-15
Curved LLT – 15 cm Tapered Core PTFE-Coated Stainless Steel with Heparin Coating	
G00575	TSCFNB-35-145-3-BH
G00573	TSCFNB-35-145-15-BH
G05624	TSCFNB-35-180-SPENCER-MOD-BH

Rosen

FIXED-CORE WIRE GUIDE

Used for angioplasty catheter exchange.

- **Heavy duty mandril** delivers the body needed for catheter exchange
- **2 cm flexible tip** for atraumatic catheter exchange
- Tightened J-shaped tip configuration allows for minimal extension beyond catheter tip
- PTFE-coated stainless steel



Order Number	Reference Part Number
Curved	
PTFE-Coated Stainless Steel	
G01774	THSCF-35-80-1.5-ROSEN
G01261	THSCF-35-145-1.5-ROSEN
G01264	THSCF-35-180-1.5-ROSEN
Curved	
PTFE-Coated Stainless Steel with Heparin Coating	
G01382	THSCF-35-145-1.5-ROSEN-BH
G01428	THSCF-35-180-1.5-ROSEN-BH


Order Number	Reference Part Number
Curved – Exchange	
PTFE-Coated Stainless Steel	
G01862	THSCF-25-260-1.5-ROSEN
G01623	THSCF-35-220-1.5-ROSEN
G01253	THSCF-35-260-1.5-ROSEN
Curved – Exchange	
PTFE-Coated Stainless Steel with Heparin Coating	
G02246	THSCF-35-260-1.5-ROSEN-BH

Double Flexible-Tipped

FIXED-CORE WIRE GUIDE

Used for access and catheterization procedures.

- Provides both **straight and 3 mm J-shaped** flexible tips on one wire
- Permits **alternative use of both ends** depending on procedural needs
- PTFE-coated stainless steel



Order Number	Reference Part Number
Straight and Curved	
PTFE-Coated Stainless Steel	
G00364	TDOC-35-50-0-3
G01306	TDOC-35-145-0-1.5
G00362	TDOC-35-145-0-3
G03794	TDOC-35-260-0-3


Order Number	Reference Part Number
Straight and Curved	
PTFE-Coated Stainless Steel with Heparin Coating	
G00363	TDOC-35-145-0-3-BH

Movable Core

WIRE GUIDE


Used to negotiate tortuous vessels.

- **PTFE-coated movable core** provides smooth advancement or withdrawal of the movable core
- Tapered distal end provides **gradual flexible-tip transition**
- Available in **straight and curved tip** configurations



Order Number	Reference Part Number
Curved – Tefcor™ Tapered	
PTFE-Coated Stainless Steel	
G01730	TCMTNA-35-145-1.5
G01678	TCMTNA-35-145-3
G01768	TCMTNA-35-145-15
G01724	TCMTNA-38-145-3
Curved – Tefcor Tapered	
PTFE-Coated Stainless Steel with Heparin Coating	
G01720	TCMTNA-35-145-1.5-BH
G01716	TCMTNA-35-145-3-BH
G01717	TCMTNA-38-145-3-BH
Straight – Tefcor Tapered	
PTFE-Coated Stainless Steel	
G01677	TMT-35-145
Straight – Tefcor Tapered	
PTFE-Coated Stainless Steel with Heparin Coating	
G01663	TMT-35-145-BH

Order Number	Reference Part Number
Amplatz Tapered	
PTFE-Coated Stainless Steel	
G03248	TCMTNA-35-145-3-AMC
Amplatz Tapered	
PTFE-Coated Stainless Steel with Heparin Coating	
G03137	TCMTNA-35-145-3-AMC-BH
Coons-Bentson	
PTFE-Coated Stainless Steel	
G27041	TSMB-35-145



Specialty Wire Guides


Whether it's to cross an infrapopliteal lesion or to drain an abscess, these specialized wire guides have been engineered to help interventionalists achieve success in specific procedures or difficult circumstances.

Approach® CTO

MICROWIRE GUIDE

Used to facilitate delivery of percutaneous catheters into the peripheral vasculature.


- Supportive **stainless steel mandril** helps in navigating support catheters smoothly through the vessel
- Wide selection of tip loads for crossing challenging lesions: **6, 12, 18 and 25 grams**
- **Shorter taper** provides enhanced pushability and penetration power



Order Number	Reference Part Number
6 gram Tip Load	
G50782	CMW-14-135-6G
G50786	CMW-14-190-6G
G50790	CMW-14-300-6G
12 gram Tip Load	
G50783	CMW-14-135-12G
G50787	CMW-14-190-12G
G50791	CMW-14-300-12G

Order Number	Reference Part Number
18 gram Tip Load	
G50784	CMW-14-135-18G
G50788	CMW-14-190-18G
G50792	CMW-14-300-18G
25 gram Tip Load	
G50785	CMW-14-135-25G
G50789	CMW-14-190-25G
G50793	CMW-14-300-25G

Tip Load Options




Approach® Hydro ST

MICROWIRE GUIDE

Used to facilitate delivery of percutaneous catheters into the peripheral vasculature.

- Supportive **stainless steel mandril** helps in navigating support catheters smoothly through the vessel
- **Shorter taper** with **hydrophilic coating** for smooth advancement through the vessel
- **Platinum coils** in the distal tip provide extra visibility



Order Number

Reference Part Number

Straight

Hydrophilic Coating


G52937	HMW-14-135-ST
G52938	HMW-14-190-ST
G52939	HMW-14-300-ST

Approach® Pro-LT

MICROWIRE GUIDE

Used to facilitate delivery of percutaneous catheters into the peripheral vasculature.

- **Linear-elastic nitinol** provides high-performance maneuverability, support and kink resistance
- **Palladium-coiled** distal tip offers extra visibility
- Available in **preformed angled** or **shapable tip** versions



Order Number

Reference Part Number

Shapable Tip

G47732	CMW-14-135
G47733	CMW-14-190
G47734	CMW-14-300
G51952	CMW-18-135
G51954	CMW-18-190
G51956	CMW-18-300

Order Number

Reference Part Number

Angled Tip

G47735	CMWA-14-135
G47736	CMWA-14-190
G47737	CMWA-14-300
G51953	CMWA-18-135
G51955	CMWA-18-190
G51957	CMWA-18-300


Not all part numbers shown on this product sheet may be approved for sale in all jurisdictions. Consult with your local Cook representative or customer service center for details.

Shorty®

INTERVENTIONAL WIRE GUIDE

Used for small vessel procedures, angioplasty and accessing distal branches or lesions in the vasculature.

- **Heat-treated stainless steel** provides 1:1 torque response
- **Shorter taper** length from mandril to tip provides distal rigidity
- Distal **platinum spring-coil tip** for increased visualization and control



Order Number	Reference Part Number
--------------	-----------------------

Straight

PTFE-Coated Stainless Steel


G11030	TPMG-16-60-SST
G11029	TPMG-16-180-SST
G11027	TPMG-18-180-SST

Coons

INTERVENTIONAL WIRE GUIDE

Used for interventional biliary drainage, abscess drainage and urology procedures.

- **Heavy-duty shaft** for support during procedures
- Gradual transition to **flexible distal tip** provides atraumatic negotiation of tortuous anatomy
- **PTFE-coated** stainless steel



Order Number	Reference Part Number
--------------	-----------------------

Straight

PTFE-Coated Stainless Steel

G12268	THSF-35-80-COONS
G02356	THSF-35-145-COONS
G02621	THSF-35-180-COONS
G02622	THSF-35-260-COONS
G04494	THSF-38-80-COONS
G02462	THSF-38-100-COONS
G02324	THSF-38-145-COONS
G02494	THSF-38-180-COONS
G02442	THSF-38-260-COONS


Not all part numbers shown on this product sheet may be approved for sale in all jurisdictions. Consult with your local Cook representative or customer service center for details.

Roadrunner® Extra-Support

MANDRIL WIRE GUIDE

Used for complex diagnostic and interventional procedures where extra support is needed for catheter exchange or manipulation.

- **Heavy-duty nitinol mandril** provides support and 1:1 torque response to distal platinum spring coil
- **Angled tip** facilitates directional control
- **PTFE-coated** mandril portion of wire guide provides low coefficient of friction for smooth manipulation



Order Number	Reference Part Number
PTFE-Coated Nitinol with Platinum Tip	
G08777	RSTF-14-180
G07557	RSTF-18-180


Order Number	Reference Part Number
Exchange PTFE-Coated Nitinol with Platinum Tip	
G08801	RSTF-14-300
G07584	RSTF-18-270
G07640	RSTF-18-300

Graduate™

MEASURING WIRE GUIDE

Used to determine accurate sizing of vessel lumen prior to angioplasty, embolization, aortic graft surgery and other interventions.

- **Gold radiopaque markers** delineate 25 cm of length for precise measuring
- **Six distal markers** are spaced 1 cm apart
- **Four proximal markers** are spaced at 5 cm increments




Order Number	Reference Part Number
Straight PTFE-Coated Stainless Steel	
G09452	MG-35-145-GRAD
G09508	MG-35-180-GRAD

McNamara

RENAL EXCHANGE WIRE GUIDE

Used for renal vascular interventional procedures.

- **Preformed curved mandril** and **short, spring-coil tip** enable wire to anchor across ostial lesions in vessels with acutely angled takeoffs
- **Nitinol mandril** provides superb kink resistance
- **PTFE coating** facilitates tracking of angioplasty catheters and other interventional devices
- **Short taper** mandril provides support to soft, **platinum spring-coil** distal tip



Order Number	Reference Part Number
--------------	-----------------------

Curved

PTFE-Coated Nitinol


G11364	TPCMG-18SP-180-MR90
G11615	TPCMG-18SP-260-MR90

Lunderquist-Ring

TORQUE INTERVENTIONAL WIRE GUIDE

Used for nonvascular interventional drainage procedures.

- **Stiff shaft** with the ability to shape and rotationally control the distal tip
- Available in **stainless steel** and **PTFE-coated** stainless steel



Order Number	Reference Part Number
--------------	-----------------------

Straight

PTFE-Coated Stainless Steel


G27031	THSF-35-125-THG
G27033	THSF-35-145-THG
G27037	THSF-38-125-THG
G01541	THSF-38-145-THG

Reuter

TIP-DEFLECTING WIRE GUIDE

Used for curving or deflecting catheter tips during selective and superselective angiography, pulmonary angiography and other applications.

- Controlled **deflection** of the wire-guide tip within the catheter lumen with strict imaging control facilitates catheter tip movement
- Distal tip of wire guide must never extend beyond distal catheter tip
- The use of a one-size-larger diameter tip-deflection wire guide will prevent wire-guide tip from inadvertently exiting catheter tip during manipulation



Order Number	Reference Part Number
Used with Tip-Deflecting Handle	
Stainless Steel	
G00763	TDW-25-100-5
G00822	TDW-25-125-5
G01162	TDW-25-145-10
G09144	TDW-25-180-5
G00384	TDW-35-100-5
G00382	TDW-35-100-10
G28700	TDW-35-120-5


Order Number	Reference Part Number
Disposable	
Stainless Steel	
G26957	DTDW-35-110-10
G07899	DTDW-35-145-10
G26958	DTDW-38-110-10
Disposable	
Stainless Steel with Heparin Coating	
G07751	DTDW-25-80-5-BH
G08159	DTDW-25-145-5-BH
G07672	DTDW-35-145-5-BH

Cope

MANDRIL WIRE GUIDE

Used for complex diagnostic and interventional procedures.

- **Soft, angled, platinum coil** distal tip for increased visualization and precise directional control
- Available with **stainless steel** or **nitinol** mandrils



Order Number	Reference Part Number
Stainless Steel	
G09593	PMG-18SP-40-COPE
G05183	PMG-18SP-60-COPE
G05822	PMG-18SP-100-COPE
G10077	PMG-18SP-125-COPE


Order Number	Reference Part Number
Nitinol – Standard Taper	
G08427	PMG-18SP-60-COPE-NT
G08687	PMG-18SP-100-COPE-NT
G08733	PMG-18SP-125-COPE-NT
Nitinol – Short Taper	
G09395	PMG-18SP-60-COPE-NT-ST

Torq-Flex®

MANDRIL WIRE GUIDE

Used for small vessel procedures, angioplasty and access of distal branches in the vasculature.


- Mandril construction provides **torque control**



Order Number	Reference Part Number
Straight Stainless Steel	
G07304	STF-18-40-AUST
G07303	STF-18-60-AUST

Order Number	Reference Part Number
Straight PTFE-Coated Stainless Steel	
G32451	THTF-18SP-150-AUST

Not all part numbers shown on this product sheet may be approved for sale in all jurisdictions. Consult with your local Cook representative or customer service center for details.



Support Wire Guides


When extra support is desired, these PTFE-coated, stainless steel wire guides offer a variety of mandril stiffness and tip configuration options.

Amplatz

SUPPORT WIRE GUIDE

Used for access and catheterization procedures.

- **Increased inner diameter** of wire guide coil provides support with a variety of mandril stiffness options
- Gradual transition to **flexible distal tip** helps facilitate atraumatic negotiation of tortuous anatomy
- **PTFE-coated** stainless steel



Order Number	Reference Part Number
Straight – Stiff	
PTFE-Coated Stainless Steel	
G27032	THSF-35-145-ASG
G28781	THSF-35-180-ASG
G27035	THSF-35-260-ASG
G27038	THSF-38-145-ASG
G27040	THSF-38-180-ASG
PTFE-Coated Stainless Steel with Heparin Coating	
G28779	THSF-35-145-ASG-BH
G27036	THSF-35-260-ASG-BH
Straight – Extra-Stiff	
PTFE-Coated Stainless Steel	
G03279	THSF-25-145-AES
G03331	THSF-25-180-AES
G04047	THSF-25-260-AES
G04003	THSF-35-80-AES
G03095	THSF-35-145-AES
G03330	THSF-35-180-AES
G03562	THSF-35-260-AES
G04001	THSF-35-300-AES
G03138	THSF-38-80-AES
G03031	THSF-38-145-AES
G03041	THSF-38-180-AES
G03579	THSF-38-300-AES
G04287	THSFNB-35-145-AES
G04367	THSFNB-38-145-AES
PTFE-Coated Stainless Steel with Heparin Coating	
G04048	THSF-35-80-AES-BH
G03290	THSF-35-180-AES-BH
G03373	THSF-35-260-AES-BH
G04187	THSF-38-80-AES-BH
G03044	THSF-38-145-AES-BH

Order Number	Reference Part Number
Curved – Extra-Stiff	
PTFE-Coated Stainless Steel	
G03584	THSCF-32-180-3-AES
G03750	THSCF-32-260-3-AES
G04237	THSCF-35-80-3-AES
G03565	THSCF-35-180-3-AES
G03460	THSCF-35-260-3-AES
G04310	THSCF-38-60-3-AES
G03540	THSCF-38-100-3-AES
G03047	THSCF-38-145-3-AES
G03217	THSCF-38-260-3-AES
PTFE-Coated Stainless Steel with Heparin Coating	
G04149	THSCF-35-80-3-AES-BH
G03547	THSCF-35-260-3-AES-BH
G03327	THSCF-38-260-3-AES-BH
Straight – Ultra-Stiff	
PTFE-Coated Stainless Steel	
G10168	THSF-35-80-AUS
G09954	THSF-35-145-AUS
G09956	THSF-35-180-AUS
G19854	THSF-35-260-AUS2
G10169	THSF-38-80-AUS
G09957	THSF-38-145-AUS
PTFE-Coated Stainless Steel with Heparin Coating	
G10051	THSF-35-145-AUS-BH
G09955	THSF-35-180-AUS-BH
G10524	THSFNA-35-145-AUS-BH


Not all part numbers shown on this product sheet may be approved for sale in all jurisdictions. Consult with your local Cook representative or customer service center for details.

Lunderquist®

SUPPORT WIRE GUIDE

Used for complex diagnostic and interventional procedures.

- Proximal portion of wire guide provides **maximum support**
- Variety of **highly flexible tip** configurations
- **PTFE-coated** stainless steel




Order Number	Reference Part Number
Straight Exchange	
PTFE-Coated Stainless Steel	
G31453	TSMG-35-260-LES
G45352	TSMG-35-300-LES
Extra-Stiff	
PTFE-Coated Stainless Steel	
G49224	TSMG-35-90-4-LES
G49225	TSMG-35-90-7-LES
G49227	TSMG-35-145-7-LES
G46729	TSMG-35-180-4-LES
G49228	TSMG-35-180-7-LES

Order Number	Reference Part Number
Curved Exchange	
PTFE-Coated Stainless Steel	
G45353	TSCMG-35-260-7-LES
G45407	TSCMG-35-300-7-LES
Double-Curved Exchange	
PTFE-Coated Stainless Steel	
G45208	TSCMG-35-260-LESDC
G45408	TSCMG-35-300-LESDC

Lunderquist®

SUPPORT WIRE GUIDE



Guide with optimal support.

SUPPORT

A PTFE-coated, stainless steel mandril is the key to Lunderquist providing maximum support and low surface friction during complex diagnostic and interventional procedures.

SAFETY

The smooth transition from the stiff mandril to the floppy gold tip and the flexibility of the tips aid vessel engagement. A radiopaque gold inner coil allows precise tip visibility and location during wire-guide manipulation.

COMPATIBILITY

Lunderquist is compatible with a wide variety of anatomies and needs. Various lengths and three tip configurations are available, including a double-curved tip design with a precurved distal tip that mimics the anatomy.

AORTIC INTERVENTION
CRITICAL CARE
ENDOSCOPY
INTERVENTIONAL RADIOLOGY
LEAD MANAGEMENT
PERIPHERAL INTERVENTION
SURGERY
UROLOGY
WOMEN'S HEALTH



Customer Service Centers

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com
www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,
E-mail: cau.custserv@cookmedical.com

