




Zilver 518[®]

VASCULAR SELF-EXPANDING STENT



Deploy a stent up to 10 mm in diameter through a 5 Fr sheath.

PROVEN RESULTS

The Zilver[®] vascular stent has been clinically proven to maintain 90% patency in the iliacs at two years.¹

LOW-PROFILE 5 FR SYSTEM

A reduced sheath size can help minimize access-site complications.²

DURABLE DESIGN

A retrospective review of iliac stenting found no fractures of the Zilver vascular stent.³



Zilver 518® RX

VASCULAR SELF-EXPANDING STENT



Zilver 518®

VASCULAR SELF-EXPANDING STENT

Zilver 518® is configured for **over-the-wire** delivery. Zilver 518RX® features **rapid-exchange** technology. Both stents are .018 inch wire guide compatible and in a 5 Fr, 125 cm delivery platform.

Order Number	Reference Part Number	Stent Diameter mm	Stent Length mm
125 cm Delivery System			
G49309	ZIVX5-18-125-6-20	6	20
G49310	ZIVX5-18-125-6-30	6	30
G49311	ZIVX5-18-125-6-40	6	40
G49312	ZIVX5-18-125-6-60	6	60
G49313	ZIVX5-18-125-6-80	6	80
G49314	ZIVX5-18-125-7-20	7	20
G49315	ZIVX5-18-125-7-30	7	30
G49316	ZIVX5-18-125-7-40	7	40
G49317	ZIVX5-18-125-7-60	7	60
G49318	ZIVX5-18-125-7-80	7	80
G49319	ZIVX5-18-125-8-20	8	20
G49320	ZIVX5-18-125-8-30	8	30
G49321	ZIVX5-18-125-8-40	8	40
G49322	ZIVX5-18-125-8-60	8	60
G49323	ZIVX5-18-125-8-80	8	80
G49324	ZIVX5-18-125-9-20	9	20
G49325	ZIVX5-18-125-9-30	9	30
G49326	ZIVX5-18-125-9-40	9	40
G49327	ZIVX5-18-125-9-60	9	60
G49328	ZIVX5-18-125-9-80	9	80
G49329	ZIVX5-18-125-10-20	10	20
G49330	ZIVX5-18-125-10-30	10	30
G49331	ZIVX5-18-125-10-40	10	40
G49332	ZIVX5-18-125-10-60	10	60
G49333	ZIVX5-18-125-10-80	10	80

Order Number	Reference Part Number	Stent Diameter mm	Stent Length mm
125 cm Delivery System			
G43769	ZIV5-18-125-6-20	6	20
G43770	ZIV5-18-125-6-30	6	30
G43771	ZIV5-18-125-6-40	6	40
G43772	ZIV5-18-125-6-60	6	60
G43773	ZIV5-18-125-6-80	6	80
G43774	ZIV5-18-125-7-20	7	20
G43775	ZIV5-18-125-7-30	7	30
G43776	ZIV5-18-125-7-40	7	40
G43777	ZIV5-18-125-7-60	7	60
G43778	ZIV5-18-125-7-80	7	80
G43779	ZIV5-18-125-8-20	8	20
G43780	ZIV5-18-125-8-30	8	30
G43781	ZIV5-18-125-8-40	8	40
G43782	ZIV5-18-125-8-60	8	60
G43783	ZIV5-18-125-8-80	8	80
G43784	ZIV5-18-125-9-20	9	20
G43785	ZIV5-18-125-9-30	9	30
G43786	ZIV5-18-125-9-40	9	40
G43787	ZIV5-18-125-9-60	9	60
G43788	ZIV5-18-125-9-80	9	80
G43789	ZIV5-18-125-10-20	10	20
G43790	ZIV5-18-125-10-30	10	30
G43791	ZIV5-18-125-10-40	10	40
G43792	ZIV5-18-125-10-60	10	60
G43793	ZIV5-18-125-10-80	10	80

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

- 1 Jaff MR, Katzen BT. Two-year clinical evaluation of the Zilver vascular stent for symptomatic iliac artery disease. *J Vasc Interv Radiol.* 2010;21(10):1489-1494.
- 2 Nasser TK, Mohler ER 3rd, Wilensky RL, et al. Peripheral vascular complications following coronary interventional procedures. *Clin Cardiol.* 1995;18(11):609-614.
- 3 Higashiyama W, Kubota Y, Sakaguchi S, et al. Prevalence, factors, and clinical impact of self-expanding stent fractures following iliac artery stenting. *J Vasc Surg.* 2009;49(3):645-652.

Customer Service Centers

AORTIC INTERVENTION
CRITICAL CARE
ENDOSCOPY
INTERVENTIONAL RADIOLOGY
LEAD MANAGEMENT
PERIPHERAL INTERVENTION
SURGERY
UROLOGY
WOMEN'S HEALTH

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com
www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335,
E-mail: orders@cookmedical.com
Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283,
E-mail: cau.custserv@cookmedical.com

