


Find the additional support that you need to exchange catheters.



Amplatz
SUPPORT WIRE GUIDES

**NOW
AVAILABLE**
1 cm
Floppy-Tip
Version

Versatile Support

Three levels of stiffness provide support options. Cook Medical recommends the extra-stiff option for the Zilver® PTX® Drug-Eluting Peripheral Stent.

Atraumatic Design

A gradual transition from the stiff shaft to the flexible tip helps you atraumatically negotiate tortuous anatomy.

Reduced Friction

The PTFE-coated stainless steel outer coil ensures smooth manipulation and eases passage.



Amplatz

SUPPORT WIRE GUIDES

Order Number	Reference Part Number	Wire Guide Diameter inch	Wire Guide Length cm	Taper Length cm	Floppy Tip Length cm
Amplatz Stiff					
Straight – PTFE-Coated Stainless Steel					
G27032	THSF-35-145-ASG	.035	145	8	3
G28781	THSF-35-180-ASG	.035	180	8	3
G27035	THSF-35-260-ASG	.035	260	8	3
G27038	THSF-38-145-ASG	.038	145	7	3
G27040	THSF-38-180-ASG	.038	180	7	3
Straight with Heparin Coating					
G28779	THSF-35-145-ASG-BH	.035	145	8	3
G27034	THSF-35-180-ASG-BH	.035	180	8	3
G27036	THSF-35-260-ASG-BH	.035	260	8	3
Amplatz Extra-Stiff					
Straight – PTFE-Coated Stainless Steel					
G03279	THSF-25-145-AES	.025	145	8	3
G03331	THSF-25-180-AES	.025	180	8	3
G04047	THSF-25-260-AES	.025	260	8	3
G04003	THSF-35-80-AES	.035	80	7	3
G04287	THSFNB-35-145-AES	.035	145	15	3
G03095	THSF-35-145-AES	.035	145	7	3
G03330	THSF-35-180-AES	.035	180	7	3
G03562	THSF-35-260-AES	.035	260	7	3
G04001	THSF-35-300-AES	.035	300	7	3
G23637	THSF-35-80-AES-ST	.035	80	7	1
G23638	THSF-35-145-AES-ST	.035	145	7	1
G23639	THSF-35-180-AES-ST	.035	180	7	1
G23640	THSF-35-260-AES-ST	.035	260	7	1
G23641	THSF-35-300-AES-ST	.035	300	7	1
G03138	THSF-38-80-AES	.038	80	9.5	3
G03539	THSF-38-100-AES	.038	100	9.5	3
G04367	THSFNB-38-145-AES	.038	145	15	3
G03031	THSF-38-145-AES	.038	145	9.5	3
G03041	THSF-38-180-AES	.038	180	9.5	3
G03579	THSF-38-300-AES	.038	300	9.5	3
Straight with Heparin Coating					
G04048	THSF-35-80-AES-BH	.035	80	7	3
G03264	THSF-35-145-AES-BH	.035	145	7	3
G03290	THSF-35-180-AES-BH	.035	180	7	3
G03373	THSF-35-260-AES-BH	.035	260	7	3
G04187	THSF-38-80-AES-BH	.038	80	9.5	3
G03044	THSF-38-145-AES-BH	.038	145	9.5	3


The new Zilver® PTX® Drug-Eluting Peripheral Stent is compatible with these sizes.*

1 cm floppy tip version now available

* For more information on wire-guide selection, see the Zilver PTX Instructions for Use.


Order Number	Reference Part Number	Wire Guide Diameter inch	Wire Guide Length cm	Taper Length cm	Floppy Tip Length cm
Amplatz Extra-Stiff					
3 mm J Curved – PTFE-Coated Stainless Steel					
G28753	THSCF-32-60-3-AES	.032	60	8	3
G03584	THSCF-32-180-3-AES	.032	180	8	3
G03750	THSCF-32-260-3-AES	.032	260	8	3
G04237	THSCF-35-80-3-AES	.035	80	7	3
G03203	THSCF-35-145-3-AES	.035	145	7	3
G03565	THSCF-35-180-3-AES	.035	180	7	3
G03460	THSCF-35-260-3-AES	.035	260	7	3
G04310	THSCF-38-60-3-AES	.038	60	9.5	3
G03540	THSCF-38-100-3-AES	.038	100	9.5	3
G03047	THSCF-38-145-3-AES	.038	145	9.5	3
G03217	THSCF-38-260-3-AES	.038	260	9.5	3
3 mm J Curved with Heparin Coating					
G04149	THSCF-35-80-3-AES-BH	.035	80	7	3
G03547	THSCF-35-260-3-AES-BH	.035	260	7	3
G12232	THSCF-35-260-3-AESW-BH	.035	260	4	2.5
G03327	THSCF-38-260-3-AES-BH	.038	260	9.5	3
Amplatz Ultra-Stiff					
Straight – PTFE-Coated Stainless Steel					
G10168	THSF-35-80-AUS	.035	80	9.5	3
G09954	THSF-35-145-AUS	.035	145	9.5	3
G09956	THSF-35-180-AUS	.035	180	9.5	3
G23642	THSF-35-80-AUS-ST	.035	80	9.5	1
G23643	THSF-35-145-AUS-ST	.035	145	9.5	1
G23644	THSF-35-180-AUS-ST	.035	180	9.5	1
G10169	THSF-38-80-AUS	.038	80	11	3
G09957	THSF-38-145-AUS	.038	145	11	3
Straight with Heparin Coating					
G10051	THSF-35-145-AUS-BH	.035	145	9.5	3
G10524	THSFNA-35-145-AUS-BH	.035	145	17	3
G09955	THSF-35-180-AUS-BH	.035	180	9.5	3
Amplatz Ultra-Stiff 2					
Straight – PTFE-Coated Stainless Steel					
G19854	THSF-35-260-AUS2	.035	260	13.5	6.5
3 mm J Curved - PTFE-Coated Stainless Steel					
G20320	THSCF-35-260-3-AUS2	.035	260	13.5	6.5

Some products or part numbers may not be available in all markets.
Contact your local Cook representative or Customer Service for details.




AMPLATZ WIRE CONSTRUCTION AND TIP TAPER DIMENSIONS

Amplatz Ultra-Stiff High Shaft Stiffness




THSF-35-XX-AUS - Standard AUS


Amplatz Extra-Stiff Medium Shaft Stiffness




THSF-35-XX-AES - Standard AES




THSF-35-XX-AUS-ST - Short Taper



THSF-35-XX-AES-ST - Short Taper




THSF-35-XX-AUS2



THSF-35-XX-AESW - "Whisker Wire"

Amplatz Stiff Standard Shaft Stiffness



THSF-35-XX-ASG - Standard Amplatz

To protect intellectual property, these diagrams are symbolic rather than specific engineering renders. Wire construction varies between product sizes. Diagrams are not to scale.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 97251996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com

