


Combine the performance you expect with increased usability.



Torque device



Scoop-shaped wire inserter



Roadrunner® UniGlide®

HYDROPHILIC WIRE GUIDE

Familiar feel

UniGlide offers a familiar feel and the performance you expect:

- The durable hydrophilic coating provides trackability and pushability.
- The tungsten-impregnated jacket enhances visibility.
- The tapered nitinol mandril offers 1:1 torque response.

Enhanced usability

Each UniGlide is packaged with the following:

- A torque device that is compatible with a broad range of wire-guide diameters.
- A unique scoop-shaped accessory that simplifies wire insertion.

320 cm length

UniGlide's more than 50 configurations include the longest peripheral hydrophilic wire on the market. The 320 cm length can facilitate radial and lower-limb procedures as well as the use of 150 cm balloon angioplasty catheters.



Roadrunner® UniGlide®

HYDROPHILIC WIRE GUIDE


Order Number	Reference Part Number	Diameter inch	Length cm
Standard Shaft Straight Tip			
G56139	HPW-18-80	.018	80
G56140	HPW-18-150	.018	150
G56141	HPW-18-180	.018	180
G30499	HPW-18-260	.018	260
G30500	HPW-18-320	.018	320
G56142	HPW-25-80	.025	80
G56143	HPW-25-150	.025	150
G56144	HPW-25-180	.025	180
G30501	HPW-25-260	.025	260
G56148	HPW-35-80	.035	80
G56149	HPW-35-150	.035	150
G56150	HPW-35-180	.035	180
G30502	HPW-35-260	.035	260
G30503	HPW-35-320	.035	320
G56154	HPW-38-80	.038	80
G56155	HPW-38-150	.038	150
G56156	HPW-38-180	.038	180
G30504	HPW-38-260	.038	260

Order Number	Reference Part Number	Diameter inch	Length cm
Standard Shaft Angled Tip			
G56162	HPWA-18-80	.018	80
G56161	HPWA-18-150	.018	150
G56160	HPWA-18-180	.018	180
G30489	HPWA-18-260	.018	260
G30490	HPWA-18-320	.018	320
G56166	HPWA-25-80	.025	80
G56167	HPWA-25-150	.025	150
G56168	HPWA-25-180	.025	180
G30491	HPWA-25-260	.025	260
G56172	HPWA-35-80	.035	80
G56173	HPWA-35-150	.035	150
G56174	HPWA-35-180	.035	180
G30492	HPWA-35-260	.035	260
G30493	HPWA-35-320	.035	320
G56178	HPWA-38-80	.038	80
G56179	HPWA-38-150	.038	150
G56180	HPWA-38-180	.038	180
G30494	HPWA-38-260	.038	260

Order Number	Reference Part Number	Diameter inch	Length cm
Stiff Shaft Straight Tip			
G56151	HPWS-35-80	.035	80
G56152	HPWS-35-150	.035	150
G56153	HPWS-35-180	.035	180
G30496	HPWS-35-260	.035	260
G30497	HPWS-35-320	.035	320
G56157	HPWS-38-80	.038	80
G56158	HPWS-38-150	.038	150
G56159	HPWS-38-180	.038	180
G30498	HPWS-38-260	.038	260

Order Number	Reference Part Number	Diameter inch	Length cm
Stiff Shaft Angled Tip			
G56175	HPWAS-35-80	.035	80
G56176	HPWAS-35-150	.035	150
G56177	HPWAS-35-180	.035	180
G30506	HPWAS-35-260	.035	260
G30507	HPWAS-35-320	.035	320
G56181	HPWAS-38-80	.038	80
G56182	HPWAS-38-150	.038	150
G56183	HPWAS-38-180	.038	180
G30508	HPWAS-38-260	.038	260

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.





Customer Service

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 972519996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: cau.custserv@cookmedical.com