

Attain both ease of movement and accessibility.




Illustration by Lisa Clark

UMI-Loc[®]

LOCKABLE UTERINE MANIPULATOR INJECTOR

Control

The slide-lock handle and grip make it easy to secure the catheter and reposition it as necessary.

Convenience


The curved shape of the injector allows for easy insertion and maneuverability.

Practicality

The single-use device is affordable and versatile.



www.cookmedical.com



UMI-Loc®

LOCKABLE UTERINE MANIPULATOR INJECTOR

Used for manipulation or repositioning of the uterus during diagnostic and surgical procedures, and for the injection of medium into the uterine cavity for diagnostic purposes.

Features and benefits:

- The slide-lock handle makes it easy to secure the injection catheter.
- The device is built for manipulative control.
- The device has an open lumen for injecting into the uterine cavity.
- The balloon prevents the device from becoming dislodged.
- Catheter markings allow for accurate placement.
- The device is designed for quick, intraoperative repositioning.
- The device contains no latex.

Cook Medical—Developing simple solutions for reproductive health.

Order Number	Reference Part Number	Diameter mm	Length cm	Balloon Volume mL
UMI-Loc				
G52974	J-UMIL-050037	5.0	37	10
UMI (Not pictured. The product is similar to UMI-Loc but lacks the slide-lock feature.)				
G52978	J-UMI-045033	4.5	33	10

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com

