

A detailed view of the Ascend AQ Ureteral Dilation Balloon Catheter. The device features a long, thin, flexible catheter with a red handle and a clear balloon at the tip. The handle has a textured grip and a clear window showing the balloon's position. The catheter is shown in a curved position, highlighting its flexibility.

Optimize dilation with a high-pressure, low-friction balloon.

AscendTM AQ[®]

URETERAL DILATION BALLOON CATHETER

High Pressure

Unique nylon balloon material has the strength to yield consistent balloon pressures up to 20 atm across the entire Ascend balloon range.

Visibility

Platinum bands provide exceptional radiopacity to ensure precise balloon dilation positioning in even the most demanding cases.

Access

AQ hydrophilic coating creates a low-friction surface for effortless advancement. The 5 Fr tapered tip aids passage through difficult anatomy while allowing placement over a .038 inch wire guide.



www.cookmedical.com



Ascend™ AQ®

URETERAL DILATION BALLOON CATHETER

The Ascend™ AQ® Ureteral Dilation Balloon Catheter combines hydrophilic coating, high visibility and high-pressure dilation to create a system designed for ease of use.

- 20 atm rated burst pressure
- AQ hydrophilic coating
- 6 Fr, 65 cm long catheter accepts a .038 inch diameter wire guide

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Balloon Length cm
Set (includes inflation device and balloon catheter)			
G32834	AUBS-4-4	4	4
G32835	AUBS-4-6	4	6
G32836	AUBS-4-10	4	10
G32837	AUBS-5-4	5	4
G32838	AUBS-5-6	5	6
G32839	AUBS-5-10	5	10
G32840	AUBS-6-4	6	4
G32841	AUBS-6-6	6	6
G32842	AUBS-6-10	6	10
G32850	AUBS-7-4	7	4
G32851	AUBS-7-6	7	6
G32845	AUBS-7-10	7	10
G32844	AUBS-10-4	10	4

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Balloon Length cm
Balloon Catheter Only			
G32825	AUB-4-4	4	4
G32826	AUB-4-6	4	6
G32827	AUB-4-10	4	10
G32828	AUB-5-4	5	4
G32829	AUB-5-6	5	6
G32830	AUB-5-10	5	10
G32831	AUB-6-4	6	4
G32832	AUB-6-6	6	6
G32833	AUB-6-10	6	10
G32847	AUB-7-4	7	4
G32848	AUB-7-6	7	6
G32849	AUB-7-10	7	10
G32843	AUB-10-4	10	4

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers - Urology

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com
www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com

