




Roadrunner[®] PC

HYDROPHILIC WIRE GUIDE




Navigate with confidence.

Navigating anatomy can be difficult. Choosing the right wire guide can also be a challenge. Cook Medical's Roadrunner® PC Wire Guide can be used to navigate a kinked or tortuous ureter for the positioning or exchange of a catheter, and can be used to bypass a large stone en route to the kidney to establish a tract. Roadrunner has a hydrophilic coating for a low-friction surface, one-to-one torque for precise control, a platinum tip for high visualization and a nitinol core for maximum deflection and kink resistance.

Roadrunner is available as either a standard-body or stiff-body wire guide. Additional options include a variety of tip configurations: angled, straight or double flexible. The double flexible tip allows for the safe introduction of the wire guide into the patient at one end while guarding against damage to scopes and other instruments placed over the opposite end.

Cook Medical—Pioneering urological products for a physician to use and a patient to trust.



Global Product Number	Order Number	Diameter inch	Length cm	Tip Configuration	Tip Length cm	Coating
Standard Body Wire Guides						
G18143	RFPC-035145	.035	145	angled	7	Slip-Coat®
G18145	RFPC-038145	.038	145	angled	7	Slip-Coat
G17540	RFSPC-035145	.035	145	angled	3	Slip-Coat
G17541	RFSPC-038145	.038	145	angled	3	Slip-Coat
G18142	RFPC-035145-0	.035	145	straight	7	Slip-Coat
G18144	RFPC-038145-0	.038	145	straight	7	Slip-Coat
G17542	RFSPC-035145-0	.035	145	straight	3	Slip-Coat
G17539	RFSPC-038145-0	.038	145	straight	3	Slip-Coat
Ink Marks						
G18824	RFSPC-035145-0-I	.035	145	straight	3	
G17898	RFSPC-038145-0-I	.038	145	straight	3	
Ink Marks and AQ® Coating						
G18629	RFSPC-035145-0-I-AQ	.035	145	straight	3	AQ
G17866	RFSPC-038145-0-I-AQ	.038	145	straight	3	AQ
Stiff Body Wire Guides						
G18154	RPC-035145	.035	145	angled	8	Slip-Coat
G18156	RPC-038145	.038	145	angled	8	Slip-Coat
G18150	RLPC-035145	.035	145	angled	16	Slip-Coat
G18152	RLPC-038145	.038	145	angled	16	Slip-Coat
G18153	RPC-035145-0	.035	145	straight	8	Slip-Coat
G18155	RPC-038145-0	.038	145	straight	8	Slip-Coat
G18149	RLPC-035145-0	.035	145	straight	16	Slip-Coat
G18151	RLPC-038145-0	.038	145	straight	16	Slip-Coat
Double Flexible with Ink Marks in 5 cm Increments						
G18791	RSPC-035145-DF-I	.035	145	double flexible		
G18793	RSPC-038145-DF-I	.038	145	double flexible		
AQ Coating						
G18792	RSPC-035145-DF-I-AQ	.035	145	double flexible		AQ
G18794	RSPC-038145-DF-I-AQ	.038	145	double flexible		AQ

Slip-Coat is a registered trademark of Angiotech Biocoatings Corp.



www.cookmedical.com

COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 USA
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll-Free Fax: 800.554.8335

COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905.640.7110, Toll Free: 800.668.0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.

95 Brandl Street, Brisbane Technology Park, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

WILLIAM COOK EUROPE ApS

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2010

URO-BM-ROADR-EN-201008

AORTIC
INTERVENTION

CRITICAL
CARE

ENDOSCOPY

INTERVENTIONAL
RADIOLOGY

LEAD
MANAGEMENT

PERIPHERAL
INTERVENTION

SURGERY

UROLOGY

WOMEN'S
HEALTH