



DAWSON-MUELLER DRAINAGE CATHETER



Now .035 wire guide compatible.

EASY TO USE

Intro-Tip™ and AQ® hydrophilic coating aid insertion and help prevent the catheter from accioning.

EASY TO SECURE

Low-profile Mac-Loc® provides firm fixation of the 10 mm pigtail.


EASY TO TRUST

Made from quality materials and inspired by clinician insight to create a premier drainage device.



DAWSON-MUELLER DRAINAGE CATHETER

ORDERING INFORMATION



Global Product Number	Order Number	Catheter Fr	Accepts Wire Guide Diameter inch	Catheter Length cm	Catheter Sideports	Trocar Stylet Needle gage
G51594	ULT6.3-35-15-P-5S-CLDM-HC	6.3	.035	15	5	19 UTW
G51595	ULT6.3-35-25-P-5S-CLDM-HC	6.3	.035	25	5	19 UTW



www.cookmedical.com

COOK MEDICAL INCORPORATED
P.O. Box 4195, Bloomington, IN 47402-4195 U.S.A.
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

COOK (CANADA) INC.
111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905.640.7110, Toll Free: 800.668.0300

WILLIAM A. COOK AUSTRALIA PTY. LTD.
95 Brandl Street, Brisbane Technology Park, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

WILLIAM COOK EUROPE ApS
Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2009 IR-BM-ULTDS-EN-200904

AORTIC
INTERVENTION

CARDIOLOGY

CRITICAL
CARE

ENDOSCOPY

INTERVENTIONAL
RADIOLOGY

PERIPHERAL
INTERVENTION

SURGERY

UROLOGY

WOMEN'S
HEALTH