




DISPOSABLE HYSTEROSCOPIC POLYP SNARE

Achieve polyp removal
with more control.



CONTROL

Duckbill shape of open snare allows enhanced control and capture of pedunculated uterine polyps.

CONVENIENCE

One-piece unit enables rapid deployment, easy capture, and quick excision in one maneuver.

PRACTICALITY

Single-use device is disposable and simple to operate.

COOK[®]
MEDICAL

DISPOSABLE HYSTEROSCOPIC POLYP SNARE



Global Product Number	Order Number	Fr	Length cm	Snare Diameter mm
G54574	J-DHPS-555015	5.5	50	15
G54573	J-DHPS-555025	5.5	50	25

The Disposable Hysteroscopic Polyp Snare (DHPS) is used to cut and coagulate polyps and fibroids for removal from the uterus under direct vision.

- Snare and handle as a single unit.
- Disposable. Intended for one-time use.
- Can be used through the working channel of a rigid or flexible hysteroscope.
- Can be used with or without electrocautery.
- Handle accepts a 2 mm monopolar plug for connection to a power source.

Cook Medical—Developing simple solutions for women’s health.



www.cookmedical.com

COOK MEDICAL INCORPORATED
P.O. Box 4195, Bloomington, IN 47402-4195 USA
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll-Free Fax: 800.554.8335

COOK IRELAND LTD.
O’Halloran Road, National Technology Park
Limerick, IRELAND, Phone: +353 61 334440, Fax: +353 61 334441

COOK (CANADA) INC.
111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905.640.7110, Toll Free: 800.668.0300, Fax: 905.640.6804

WILLIAM A. COOK AUSTRALIA PTY. LTD.
95 Brandl Street, Brisbane Technology Park, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

© COOK 2010 WH-BM-DHPS-EN-201006