

Manipulate and remove uterine polyps under direct vision.



NCircle[®]
NITINOL HYSTEROSCOPIC POLYP SNARE



www.cookmedical.com



NCircle®

NITINOL HYSTEROSCOPIC POLYP SNARE

The NCircle Polyp Snare provides a simple way to manipulate and remove uterine polyps under direct vision. Designed for use on smaller polyps, it is easy to operate and slides open and closed with little effort.

- The nitinol memory wire opens easily in the uterine cavity to allow the capture and removal of polyps.
- The small 5.0 Fr outer sheath fits through the working channel of many hysteroscopes.
- The snare is intended for one-time use.

Cook Medical—Developing simple solutions for women's health.

Order Number	Reference Part Number	Catheter Fr	Length cm
G18738	J-NTSE-050035-UDH-HPS	5.0	35

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 97251996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: cau.custserv@cookmedical.com

