

Maximize control
with Balloon-Assisted
Tracheostomy (BAT).



Dolphin BT™

CIAGLIA BALLOON-ASSISTED TRACHEOSTOMY INTRODUCER



www.cookmedical.com

DOLPHIN BT™ TRACHEOSTOMY SET

The Dolphin BT tracheostomy introducer is the first device to combine stoma formation and tracheal tube insertion into a simple and efficient one-step procedure: balloon-assisted tracheostomy (BAT). The Dolphin BT tracheostomy introducer is designed specifically to maximize operator control and prevent trauma to the anterior tracheal wall and surrounding structures.



The Dolphin BT™ Tracheostomy Set includes the following components:

- Dolphin BT tracheostomy introducer
- Cook Sphere® Inflation Device
- Wire guide with positioning marks
- 14.0 Fr dilator
- Disposable #15 safety scalpel
- 5" stainless steel curved hemostat
- 18 gage introducer EchoTip® needle
- 18 gage disposable PTFE catheter introducer needle

- Disposable syringe
- 36" x 48" drape with a clear plastic window
- Lubricated jelly
- Gauze sponges
- Needle holder cup
- Tyvek tape measure
- CSR wrap 30" x 30"

The set is also available in a version with a sterile tracheostomy tube.

Order Number	Reference Part Number	Loading Dilator Fr	Tracheostomy Tube ISO size	Tracheostomy Tube Jackson size
G36101	C-PTBS-2100	21.0	-	-
G36102	C-PTBS-2400	24.0	-	-
G36103	C-PTBS-2600	26.0	-	-
G36104	C-PTBS-2700	27.0	-	-
G36105	C-PTBS-2800	28.0	-	-
G36123	C-PTBS-3000	30.0	-	-

With VersaTube™ Tapered Tracheostomy Tube

G23548	C-PTBS-2100-VT7	21.0	7	-
G55593	C-PTBS-2400-VT8	24.0	8	-
G55595	C-PTBS-2700-VT9	27.0	9	-

With Shiley® Percutaneous Cuffed Tracheostomy Tube

G36116	C-PTBS-2600-PERC6	26.0	-	6
G36117	C-PTBS-2800-PERC8	28.0	-	8

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

DOLPHIN BT™ TRACHEOSTOMY TRAY



The Dolphin BT Tracheostomy Tray includes the following components:

- Dolphin BT tracheostomy introducer
- Cook Sphere® Inflation Device introducer needle
- Wire guide with positioning marks
- 14.0 Fr dilator
- Disposable #15 safety scalpel
- Curved mosquito clamps
- 5" stainless steel curved hemostat
- 18 gage introducer EchoTip® needle
- 18 gage disposable PTFE catheter
- 22 gage disposable infiltration needle
- 25 gage disposable infiltration needle
- Tinted ChloraPrep® One-Step (2)
- 5 mL vial of lidocaine epinephrine 1.5% (2)
- Double-swivel connector
- Suture with needle
- Disposable syringes
- 36" x 48" drape with a clear plastic window
- Lubricated jelly
- Gauze sponges
- Filter straw
- Tracheostomy tube holder
- Drain sponges
- Needle holder cup
- Tyvek tape measure
- CSR wrap 30" x 30"

The tray is also available in a version with a sterile tracheostomy tube.

Order Number	Reference Part Number	Loading Dilator Fr	Tracheostomy Tube ISO size	Tracheostomy Tube Jackson size
G36106	C-PTBSY-2100	21.0	-	-
G36107	C-PTBSY-2400	24.0	-	-
G36108	C-PTBSY-2600	26.0	-	-
G36109	C-PTBSY-2700	27.0	-	-
G36110	C-PTBSY-2800	28.0	-	-
G49506	C-PTBSY-3000	30.0	-	-

With VersaTube™ Tapered Tracheostomy Tube

G23551	C-PTBSY-2100-VT7	21.0	7	-
G55594	C-PTBSY-2400-VT8	24.0	8	-
G55596	C-PTBSY-2700-VT9	27.0	9	-

With Shiley® Percutaneous Cuffed Tracheostomy Tube

G36118	C-PTBSY-2600-PERC6	26.0	-	6
G49507	C-PTBSY-2800-PERC8	28.0	-	8

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Customer Service

EMEA: EDI - www.cookmedical.com/edi.do

Distributors: +353 61239240, ssc.distributors@cookmedical.com

Austria: +43 179567121, oe.orders@cookmedical.com

Belgium: +32 27001633, be.orders@cookmedical.com

Denmark: +45 38487607, da.orders@cookmedical.com

Finland: +358 972519996, fi.orders@cookmedical.com

France: +33 171230269, fr.orders@cookmedical.com

Germany: +49 6950072804, de.orders@cookmedical.com

Hungary: +36 17779199, hu.orders@cookmedical.com

Ireland: +353 61239252, ie.orders@cookmedical.com

Italy: +39 0269682853, it.orders@cookmedical.com

Netherlands: +31 202013367, nl.orders@cookmedical.com

Norway: +47 23162968, no.orders@cookmedical.com

Spain: +34 912702691, es.orders@cookmedical.com

Sweden: +46 858769468, se.orders@cookmedical.com

Switzerland - French: +41 448009609, fr.orders@cookmedical.com

Switzerland - Italian: +41 448009609, it.orders@cookmedical.com

Switzerland - German: +41 448009609, de.orders@cookmedical.com

United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283

E-mail: cau.custserv@cookmedical.com